


enFlow™

IV fluid and blood warmer

System Manual
44000024 Rev AC



About this manual

This manual has been developed to provide the user with the information necessary to operate and maintain the enFlow IV fluid/blood warming system. It is important that all medical personnel that operate this device read and understand all the information contained within this System Manual. This material is not meant as a substitute for formal training in the use of intravenous delivery systems, which may be required by local, regional or state protocol. As with any medical device, please consult your local medical director or governing agency for further information and requirements. If you have questions or concerns regarding this manual or product, please contact Customer Service or Technical Support for assistance:

Customer Service (United States)

Phone: 1.833.327.3284

For customers outside the United States, please contact your local Vyaire representative or local vendor.

Australian Sponsor:
Vyaire Medical Pty Ltd
Suite 5.03, Building C,
11 Talavera Road,
Macquarie Park, NSW 2113



MANUFACTURED FOR
Vyaire Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045 USA
1-833-327-3284
customersupport@vyaire.com
www.vyaire.com



AUTHORIZED
REPRESENTATIVE IN THE
EUROPEAN COMMUNITY
EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands











































2797

© 2020 Vyaire. Vyaire, the Vyaire logo and all other trademarks or registered trademarks are property of Vyaire Medical, Inc., or one of its affiliates.

Symbols used on the equipment

The following symbols may be viewed on any of the products or accessories that comprise the enFlow IV fluid/blood warming system.

Symbol	Symbol description	Symbol	Symbol description
	Batch Code		On
	Catalog Number		Off
	Serial Number		Not Made With Natural Rubber Latex
	National Stock No. (US Military)		Expiration Date
	Single Use Only; Do Not Re-Use		Direct Current
	Sterilized Using Irradiation		Alternating Current
	Keep Dry		Type BF Applied Part, Defibrillation-Proof
	Do Not Re-Sterilize		Do Not Use if Package is Damaged
	Caution		Fuse
	Temperature; Thermometer		Non-Pyrogenic
	Danger High Voltage	Note 	This symbol indicates that additional information is being provided.

Symbol	Symbol description	Symbol	Symbol description
	Electric Energy		Effect or action in both directions away from reference point (Open)
	Temperature Range		Effect or action in both directions towards a reference point (Close)
	Not made with Di(2-ethylhexyl) phthalate		In transport applications it is advised to cushion and insulate the Warmer from the patient's skin and apply the Warmer as loosely as acceptable checking regularly for signs of potential pressure-related injury
	Manufacturer	IP67	Degree of protection provided by enclosure, dust tight, temporary water immersion
IP31	Degree of protection provided by enclosure, no ingress of object > 2.5 mm diameter, protected against dripping water	IP68	Degree of protection provided by enclosure, dust tight, continuous water immersion
	Caution: U.S. Federal law restricts this device to sale by or on the order of a physician		Surface may be hot
	Consult Instructions for Use		Do not encase the Warmer with any external coverings like: towels, sheets, blankets or drapes. Covering the Warmer restricts the natural convection of heat.
	System Fault XX		Low Battery
	Lock or Password Required		Unlock
	Do Not Throw in Trash	CE	The CE Mark is the manufacturer's or importer's mark of conformity declaring compliance with all applicable directives (Safety, EMC, Machinery, Medical and others).
	System Fault XX		Press Any Key to Mute Alarm








Symbol	Symbol description	Symbol	Symbol description
	Mute the Audible 'High Priority Alarm' For 1 Minute		Not Heating
	Medical device		Humidity limitation
	Atmospheric Pressure		Authorised Representative in the European Community
	Intertek is accredited by OSHA as a NRTL, as well as by the Standards Council of Canada. This mark indicates that the product has been tested to CAN/CSA C22.2#60601-1:2014 Ed 3, AAMI ES60601-1:2005 + A1:2012, IEC 60601-1:2005 Ed 3 + A1:2012, IEC 60601-1-6:2010 Ed 3 + A1:2013, IEC 60601-1-8:2006 Ed 2 + A1:2012, IEC 60601-1-9:2009 Ed 1 + A1:2013, IEC 62304:2006 Ed 1 + A1:2015, IEC 62366:2007 Ed 1 + A1:2012		

Table of contents

Warnings	7
Cautions	7
enFlow IV fluid/blood warming system description	8
Indications for use	9
Clinical and training information	9
Unpacking the enFlow IV fluid/blood warming system	9
To begin operation of the enFlow IV fluid/blood warming system	10
enFlow Controller (Model 121 series) indicators and operation	12
Controller (Model 121 series)—Setup Instructions	12
Controller Display	13
enFlow Warmer (Model 100 series) indicators and operation	13
Cleaning the enFlow IV fluid/blood warming system components	15
Caution	15
Warmer and Controller Cleaning Procedure	15
Warmer High Level Disinfection Procedure	15
Warmer and Controller Intermediate Disinfection Procedure	15
Storing the enFlow IV fluid/blood warming system components	16
Servicing the enFlow IV fluid/blood warming system components	16
Instructions for replacing the Controller clock battery	16
Instructions for changing the Controller fuse	18
enFlow fluid warming system temperature controls and alarms	18
Temperature control	18
Audible/Visual alarm	19
enFlow troubleshooting	19
Electromagnetic interference	19
Interference confirmation	19
Interference reduction steps	19
enCheck™ (Model 400) user guide	20
Cleaning the enCheck	21
Appendix A: Technical specifications	22
Appendix B: Fault code table	26
Appendix C: Warming system response by temperature	27
Appendix D: Parts list	28
Appendix E: Preventative maintenance procedure	29
Appendix F: enFlow IV fluid/blood warming system operational checklist	38
Appendix G: enFlow IV fluid/blood warming system operational checklist—enCheck Model 400 method	39
Appendix H: Glossary	40

Warnings

- All IV fluid bags must be vented of air per IV fluid manufacturers' directions prior to connecting to the infusion set. Standard IV line protocols for priming the complete infusion set, the enFlow Disposable Cartridge, and the extension set must be followed before connecting to a patient. Care must be taken to ensure there is not sufficient air in the fluid bag and lines to cause an air embolism.
- The 'High Priority Alarm' is a flashing RED LED, a flashing RED controller display, and an audible alarm, indicating that the infusate is over temperature. Stop the fluid flow, and slide the Warmer covers open to stop warming. If the above occurs, then replace the Warmer and contact Technical Support. The attending practitioner should remain within 4 m of the patient when the device is in use to enable visualization of the enFlow display and hear the audible high priority alarm.
- The Warmer contains magnets; do not operate within 15 cm (6 in) of a pacemaker or other devices that may be sensitive to strong magnetic fields.
- The Disposable Cartridge may be a potential biohazard during or after use. Handle and dispose of in accordance with acceptable medical practice and applicable regulations.
- Do not use in the presence of flammable anesthetics.
- Replace the fuses with Bussmann® part #S500-5-R or equivalent.
- The Disposable Cartridge should not be used for more than 24 hours.
- Ensure that the Disposable Cartridge expiration date has not passed.
- If the IV line runs dry, disconnect the Disposable Cartridge from the Warmer. Re-prime the entire IV system using aseptic techniques. Ensure all the air is removed from both the line and the Disposable Cartridge. Replace the Disposable Cartridge in the Warmer.
- The enFlow Warmer is to be used only with approved enFlow power sources and the enFlow Disposable Cartridge.
- To avoid risk of electrical shock, this equipment must only be connected to a supply main that is grounded. Should the need arise the device may be disconnected by appliance coupler.
- Do not connect the enCheck to an enFlow system while it is in any way connected to patient.
- Do not Service or perform maintenance while in use with the patient. Before performing maintenance, disconnect the enflow system from the patient, shut it down, and disconnect it from AC power. All operator maintenance must be performed with the patient off the enflow system. Failure to do so can result in electric shock to the patient and operator.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the enFlow controller and warmer, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Cautions


- Please contact your sales rep or customer service to make sure you are using the latest revision designs.
- Follow the AABB "Guidelines for the Use of Blood Warming Devices" (© 2002) which caution against warming when administering platelets, cryoprecipitate, or granulocyte suspensions.
- Some drugs or drug preparations may be sensitive to warming. As with any fluid or blood warming system, carefully review the drug manufacturer's literature for information about thermal sensitivity.
- The disposable cartridge contains polycarbonate and Parylene in the fluid path. Review the preparation or solution manufacturer's instructions for use about chemical sensitivity.
- Do not affix, place or bind the Warmer directly to a patient during general use.
- Do not wrap the Warmer in towels, sheets, blankets or drapes.
- Insulated Warmer Strap (980304VS30) may be used with the Warmer for transport purposes only, such as field use (for the military) and inter-hospital transport and ambulatory circumstances. Do not strap for non-transport scenarios. If the enflow system is used for pre-hospital transport or transfer to another facility and there is a desire to secure the warmer with more than the cord clip then the clinician must implement the following instructions:
 - ▶ Stabilize the warmer using the 980304VS30 insulated warmer strap, do not use any other strap or towels to secure the warmer to the patient.
 - ▶ Place an insulating and cushioning fabric layer, such as soft cotton towels or gauze, at least .25" or 6 mm thick in between the underside of the warmer and the patient. Do not use foam or gel pads. Cushioning the patient from the warmer will help prevent perioperative peripheral neuropathies.
 - ▶ Attach the insulated warmer strap as loosely as possible, taking care to reduce the possibility of exsanguinating the appendage or the area of attachment. Check regularly for signs of potential pressure related injury.
- The Warmer heating surface and Disposable Cartridge can get quite warm when heating cold IV fluids/blood at high flow rates. Wait a few seconds after stopping the IV fluid/blood flow before removing the Disposable Cartridge.
- The Controller should only be plugged into a hospital grade outlet.
- Do not block the fan in the Controller as this may cause overheating.
- Although the Warmer has been tested to insure it will survive a drop of 1 m (3.28 ft), care should be taken that the device is not dropped to reduce the potential of damage.
- Do not clean with:
 - ▶ ketones (MEK, acetone, etc.)

- ▶ abrasive cleaners.
- Do not sterilize the Warmer with:
 - ▶ steam sterilization (autoclave)
 - ▶ dry heat.
- Do not sterilize the Controller.
- Do not spray or pour cleaning solutions directly on the Controller.
- Do not allow cleaning solutions to accumulate on the Controller.
- When using the Controller mounted to an IV pole, it must be tightly secured on the pole no higher than 122 cm (48 in) from the ground. The pole should have a base diameter of no less than 61 cm (24 in). A Controller mounted too high on the IV pole may cause instability. IV pole accessories or the attachment of fluid bags may also cause instability.
- Normal wear and tear during use of the Warmer may cause the device to be susceptible to fluid ingress. Carefully inspect the heating surface of the Warmer for tears or foreign matter before each use and take out of service if necessary.
- Always secure the infusion set with the provided IV Line Clip on the Warmer power cable to prevent kinking in the line.
- Do not use a stiff bristle brush or sharp probe to remove foreign material.
- Do not use compressed air to dry.
- Avoid puncturing the heating surface. If damaged, remove the Warmer from service and replace immediately.
- This equipment is not intended for use in an oxygen rich environment (defined as $>22\% \text{ O}_2$).
- No modification of this equipment is allowed.
- Do not position the device in a way that makes it difficult to disconnect the device.
- Due to highly stable components, microprocessor control, and built-in self-tests, an annual performance check is sufficient.
- The steps listed in the enFlow IV fluid/blood warming system operational checklist (Appendix F) should be performed at least once a year, or as required by your accrediting body.
- The enCheck is to be used only with the enFlow Warmer and Controller.
- Do not touch the enCheck contact plate surface during or immediately after use since it may be very hot.
- The warmer and controller may be a potential biohazard during or after use. Handle and dispose of in accordance with acceptable medical practice and applicable regulations at the end of their life.
- Dispose of the product, its accessories and internal and external battery according to all applicable local regulations after the end of expected service life.
- Waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Contact a Vyair authorized representative for information concerning the decommissioning of your equipment.


enFlow IV fluid/blood warming system description

The enFlow IV fluid/blood warming system consists of the enFlow Warmer (Model 100 series), the enFlow Controller (Model 121 series), the enFlow Disposable Cartridge with IV extension set (Model 202) or without IV extension set (Model 200). Within seconds, this warming system delivers normothermic infusate to the patient at flow rates of Keep Vein Open (KVO defined as 2 mL/min) to 200 mL/min when input fluid temperature is 20 °C.


The Warmer is the reusable heating unit designed to work in conjunction with the Disposable Cartridge. Two multicolored light emitting diode (LED) indicators on the Warmer indicate its power status and the fluid/blood infusate temperature. The infusate within the Disposable Cartridge is warmed when in contact with the heating surface of the Warmer. This surface is heated by means of electrical resistance. The Warmer contains redundant temperature sensors to help ensure fluid temperature accuracy and reliability. It also includes two independent over-heating protectors. Continuous internal diagnostics monitor essential components and system parameters when heating fluid/blood.



The Controller serves as a power supply for the Warmer unit. The Controller is designed to mount on an IV pole or sit on a table top. The front panel includes a Controller reading in degrees Celsius, as well as a keypad, which controls the clock and the mute feature. The Controller display is always shown right-side-up.



Each Disposable Cartridge and the Disposable Cartridge with IV extension set are radiation sterilized and non-pyrogenic, not made with natural rubber latex or DEHP. The Disposable Cartridge connects to the IV extension set or any infusion set employing standard luer connectors. Once primed, the Disposable Cartridge in conjunction with the Warmer and the Controller combine to complete the enFlow IV fluid/blood warming system.



Indications for use

The enFlow IV fluid/blood warming system's intended use is for warming blood, blood products and intravenous solutions prior to administration. It is designed to be used by healthcare professionals in hospital, clinical and field environments to help prevent hypothermia.


Clinical and training information

Operators must be trained to set up and deliver blood/ IV solutions in a medically approved manner, including aseptic techniques and standard hospital procedures. Use of the enFlow IV fluid/blood warming system, when properly administered, will help to prevent hypothermia and the complications arising therefrom.

Unpacking the enFlow IV fluid/blood warming system


Upon receipt of the enFlow IV fluid/blood warming system components, visually inspect the shipping containers and internal contents for damage that may have occurred during shipment. If there is any visible or mechanical damage to the contents, or if the order is incomplete, please contact Customer Service immediately. The components for each model are reflected below:

- Model 100
 - ▶ Warmer
 - ▶ Warmer Cord Clip
 - ▶ IV Line Clip
 - ▶ USB Manual
 - ▶ Warmer Mount (Warmer Mount Instructions)
 - ▶ Patient Leakage Report
 - ▶ Certificate of Conformance
- Model 121
 - ▶ Controller
 - ▶ Patient Leakage Report
 - ▶ USB Manual
 - ▶ Warmer Mount (Warmer Mount Instructions)
- Model 200
 - ▶ Disposable Cartridge
 - ▶ IFU




Model 202

- ▶ Disposable Cartridge with IV extension set
- ▶ IFU




To begin operation of the enFlow IV fluid/blood warming system


- a. Place the Controller on a firm, level surface or on an IV pole with an outside dimension of no greater than 3.0 cm (1.25 in). Two Controllers may also be mounted next to each other on an IV pole as shown below. The Controller's display will have a right-side-up orientation regardless of its position.



- b. Plug the Controller into a hospital grade outlet.
- c. Setting the clock to the local time is optional, but usually done on initial use. No changes in performance are affected by the clock's setting. (Please refer to the "enFlow Controller (Model 121 series) indicators and operations" section for directions to set the clock).
- d. Connect the Warmer cable to the Controller. This action is accomplished in three steps:
1. Insert the male plug end of the Warmer into the female receptacle on the front face of the Controller. Push it in so that the plug cover is tight against the receptacle.




2. The plug and receptacle are keyed in both orientation and configuration. This feature ensures that the Warmer can only be plugged in properly. Additionally, it prevents other plug devices from fitting into this receptacle.
 3. Turn right to lock. (See arrows on label.)
- e. The rear mounted I/O (ON/OFF) switch on the Controller turns the power on and off. Switch the Controller to ON. Upon startup, the Controller conducts a self-test. The power indicator illuminates green, the Controller display flashes "enFlow," a short audible beep occurs, and the LED's light up for about one (1) second. **Note** ⚠ : The Controller automatically switches for operation at either 100, 115, or 240 VAC (100-240 VAC).




f. Open the covers on the Warmer by sliding them apart.


g. Connect the infusion set and/or extension set to the Disposable Cartridge; then prime with fluid using standard medically approved protocols. Next, connect the infusion set to the patient and place the Disposable Cartridge into the Warmer.




h. Completely close the covers on the Warmer by pressing down on the cartridge and sliding the covers inward toward each other until they meet. Upon closing the covers, a short audible beep occurs indicating that the Warmer self-test is being performed and confirms operation of temperature sensors and alarm indicators. After this process is complete, regulated power is delivered to the Warmer's heating surface, which then begins heating the infusate through the Disposable Cartridge. Adjust the fluid flow to the desired rate.



i. Place the IV line in the IV Line Clip in order to prevent it from kinking.



j. The Warmer is designed to be placed on the bed and/or attached to patient coverings in close proximity to the site of infusion using the cord clip P/N 980309VS-20. Cushion the patient from the Warmer to aid in the prevention of perioperative peripheral neuropathies or heat-related dermal injury.




k. Do not wrap the Warmer in towels, sheets, blankets or drapes.




l. Opening the Warmer covers immediately stops the **heating** but **not** the **flow**.

m. To turn off the device, use the switch located at the rear of the Controller.



enFlow Controller (Model 121 series) indicators and operation



Controller (Model 121 series)—setup instructions

1. Set the clock

- To modify the initial default mode of the Controller, press the center button of the front keys of the display window prior to pushing the power switch that is located in the back of the unit. The key must be held down continuously until the clock is displayed.
- Once the clock appears, toggle the center button to move through the fields. To change a field, press the right or left arrow keys.

2. Set the Controller display default mode to symbols (International Mode)

- Continue to toggle through to the set zone field. Press either the right or left arrow key to access the padlock symbol to enter the password screen.
- The password is 781.


- c. An underscore will display under the first digit field. Press the arrow key to set the appropriate digit.
- d. Next, press the center key again to toggle to the next place.
- e. Repeat steps three and four for the second and third digits. If an incorrect password is entered, the system reverts to the set zone screen.
- f. After the third digit is set, press the center key again. First a green padlock and then the zone USA will appear. Press the right arrow to change to INTL.
- g. Press the center key again to set. The system begins to operate.
- h. After the initial setup, whenever the Controller is powered on, the display screen defaults to the last mode entered.

Controller display

The Controller display continuously reflects the specific infusate temperature that the Warmer monitors and maintains. The various readouts that may be depicted on the Controller display are described in the following tables:

Table 1—Controller Display: Normal Operating Mode		
Activity	Display reads	Display color and function
Warmer is connected and power is engaged.	Temperature and Clock 40 °C 9:00 A	Identical to Warmer temperature LED
Warmer is not connected, but Controller is powered on.	Not Heating	Yellow
Warmer is connected, but covers are open on Warmer.	Not Heating	Yellow
Warmer is connected, and covers are either open or closed on Warmer; however, Cartridge is not in Warmer.	Not Heating	Yellow

Table 2—Controller Display: Alarm Mode		
Activity	Display reads	Display color and function
Warmer over temperature	Display alternates between Over Temp and Press Key to Mute	Identical to Warmer temperature LED
Mute button activated	Over Temp Muted	Identical to Warmer temperature LED
Fault detected	System Fault XX If a system fault message is on the Controller display, refer to the system fault section Appendix B or contact Technical Support.	Red High priority alarm

Table 3—Controller Display: Setup Mode		
Activity	Display reads	Display color and function
While powering up the Controller, hold the center button down until the clock screen is displayed. Then, use the buttons to set the clock.	09:00 A 	Blue

enFlow Warmer (Model 100 series) indicators and operation

The Warmer monitors and maintains the infusate temperature at 40 °C ± 2 °C. On the top of the Warmer, there are two status indicator lights (multicolored LEDs), which reflect the following:

Power - Indicates the power and operational status of the Warmer.

Temperature - Indicates that the infusate temperature is within an acceptable operating range (35 °C to 42 °C).




Table 4: Warmer indicator status

Status	Warmer covers	Disposable cartridge	Power LED indicator	Temperature LED indicator	Audible indicator	Description	Action required
Ready	Open or closed	None	Flashing Green every 3 seconds	Unlit	None	Warmer unit has power, but is not heating	None
Device power up	Closed	In place	Red ½ second duration	Red ½ second duration	½ second beep	Successful device power up and over temperature circuit test	Observe. If the LED does not flash red, replace the Warmer; and contact Technical Support.
In operation	Closed	In place	Solid Green	Flashing Blue	None	Infusate temperature is < 33 °C	None
In operation	Closed	In place	Solid Green	Solid Blue	None	Infusate temperature is ≥ 33 °C and < 35 °C	None
In operation	Closed	In place	Solid Green	Solid or Flashing Blue >30 seconds	None	Warmer is unable to heat the infusate within operational range. Infusate temperature is < 20 °C and/or the flow rate is > 200 mL/min.	Reduce the flow rate if possible. If there is no change in operational temperature, consider replacing the Warmer and contact Technical Support.
In operation	Closed	In place	Solid Green	Solid Green	None	Infusate temperature is ≥ 35 °C and ≤ 42 °C	None
In operation	Closed	In place	Solid Green	Solid Yellow	None	Infusate (and/or ambient temperature) is > 42 °C but less than an Over Temp condition	Observe. This state whereby the infusate is > 42 °C should only be entered periodically during changes in flow rate or infusate temperature.
In operation	Closed	In place	Flashing Red	N/A-Varies	Continuous audible bursts	Internal failure in the Warmer	Replace the Warmer if this occurs, and contact Technical Support.
Continuous operation	Closed	In place	Solid Green	Flashing Red High priority alarm	Continuous audible bursts	Infusate (and/or ambient temperature) is > 45 °C signifying an Over Temp condition	Stop the fluid flow, and slide the Warmer covers open to stop warming. Replace the Warmer if this occurs, and contact Technical Support.

Refer to Appendix C for the chart on “Warming system response by temperature.”

Refer to Warnings for additional information.

Cleaning the enFlow IV fluid/blood warming system components

Caution

- Do not clean with:
 - ▶ Ketones (MEK, acetone, etc.)
 - ▶ Abrasive cleaners
- Do not sterilize the Warmer with:
 - ▶ Steam sterilization (autoclave)
 - ▶ Dry heat
- Do not sterilize the Controller.
- Do not spray or pour cleaning solutions directly on the Controller.
- Do not allow cleaning solutions to accumulate on the Controller.

The Warmer and Controller are chemically resistant to most common hospital grade instrument cleaning solutions and non-caustic detergents. The following list of approved cleaning solutions may be used to clean the Warmer and Controller:

- 1.3 mL of Dawn® Ultra Dishwashing Liquid diluted in 1 liter of 80 – 90 ° F water
- CaviWipes®
- 4% Hydrogen Peroxide
- CIDEX OPA

Warmer and Controller Cleaning Procedure

1. Disconnect the fluid warmer from the controller.
2. Prepare a mild dishwashing cleaning solution or 4% hydrogen peroxide solution
3. Using Solution, dampen a low linting or lint-free soft-cloth. Wipe the surfaces of the Controller unit to remove debris.
Note: Do not spray or pour cleaning solutions directly on to the Controller
4. Pay particular attention to the view screen and soft-buttons on the front of the unit.
5. Repeat this step replacing the saturated soft-cloth with a fresh one if there is a soiling of the soft-cloth.
6. Repeat this step replacing with a fresh saturated soft-cloth.
7. Using water, dampen a lint-free soft cloth. Wipe the entire surfaces of the Controller unit for 30 seconds to remove any chemical residues.
8. Fluid Warmer will be completely immersed in mild dishwashing cleaning solution or 4% hydrogen peroxide solution for three (3) to five (5) minutes.
Note: Do not immerse the warmer's electrical plug connector
9. A soft bristle brush may be used to access small cracks and crevices until visibly clean.
10. Using cleaning solution, dampen a lint-free soft cloth. Wipe the surfaces of the warmer's electrical plug connector.
11. Rinse the fluid warmer thoroughly with distilled water to remove all detergent residues for one (1) minute.
12. Allow devices to air dry

Warmer High Level Disinfection Procedure:

1. Disconnect the fluid warmer from the controller.
2. Fluid Warmer will be completely immersed in the CIDEX OPA solution, for a minimum of twelve (12) minutes. The minimum exposure temperature is 20C.
Note: Do not immerse the warmer's electrical plug connector
3. Using CIDEX OPA, dampen a lint-free soft cloth. Wipe the surfaces of the electrical plug connector.
4. Following removal from CIDEX OPA solution, thoroughly rinse by immersing completely in a large volume of water (two (2) gallons). Keep the device totally immersed for a minimum of one (1) minute.
5. Repeat the rinsing procedure two (2) additional times, for a total of three (3) rinses.
6. Allow devices to air dry out of exposure to direct sunlight.

Warmer and Controller Intermediate Disinfection Procedure:

Warmer:

1. Fluid Warmer will be completely wiped with Caviwipes for three (3) to five (5) minutes.
Note: Do not immerse the warmer's electrical plug connector.
2. A soft bristle brush may be used to access small cracks and crevices until visibly clean.
3. Using CaviWipes, wipe the surfaces of the warmer's electrical plug connector.
4. Rinse the fluid warmer thoroughly with distilled water to remove all detergent residues for one (1) minute.

5. Allow devices to air dry.

Controller:

1. Disconnect the fluid warmer from the controller.
2. Using CaviWipes, wipe the surfaces of the Controller unit to remove debris.
While applying the disinfectant ensure that all surfaces of the device remain wet for at least three (3) minutes. Reapply using a fresh Caviwipe as needed.
3. Wipe the external surfaces of the Enflow Controller Unit with a fresh Caviwipe, make sure disinfectant is applied to device seams.
Note: Do not spray or pour cleaning solutions directly on to the Controller
4. Repeat this step replacing the CaviWipe with a fresh one if there is soiling of the CaviWipe.
5. Repeat this step replacing with a fresh CaviWipe.
6. Wipe the surfaces of the electrical plug connector.
7. Using water, dampen a lint-free soft-cloth. Wipe the entire surfaces of the Controller unit for 30 seconds to remove any chemical residues.

Storing the enFlow IV fluid/blood warming system components

The Warmer and Controller should be stored in a clean, dust free environment.

Servicing the enFlow IV fluid/blood warming system components

The enFlow IV fluid/blood warming system components have been designed to be durable and long lasting. The system uses current Surface Mount Technology (SMT) and materials. If a System Fault error occurs, then remove the device from service and refer to the table in Appendix B for proper evaluation and handling. Do not return the device or place it back into service unless it has been cleaned per the instructions above and evaluated by a trained technician using the instructions in Appendix E and Appendix F. If the unit stops working properly, contact Customer Service to obtain an Return Goods Authorization (RGA) number prior to returning the unit to Vyair. If damage has occurred to the heating surface, immediately remove it from service.

RGA number (United States)

The Technical Support Representative will troubleshoot your product issue with you on the phone. If it is necessary to return a product under warranty, a replacement unit will be shipped to you within 48 hours. (If the product is no longer under warranty, the Customer Advocacy Representative will discuss replacement options.) You will be issued a RGA number. You will be instructed to return the product in sufficient packaging to prevent damage in transit, clearly marking the RGA number on the outside of the box. The return address will be provided to you.

Warmer (Model 100 series)

The Warmer is permanently sealed against fluid ingress and has no user serviceable parts inside.

Controller (Model 121 series)

See "Servicing the enFlow IV fluid/blood warming system components" below for information on user serviceable parts inside the Controller. Check the fuses located in the outside power entry module if the Controller fails to function. The power cord must be removed to do this.

Replaceable parts for the Controller:

- Clock battery (See below for replacement instructions)
- Fuse – Bussmann #S500-5-R or equivalent (See below for replacement instructions)
- Pole clamp screw
- Power cord
- Warmer mount

The Controller should be subjected to routine safety checks as required by local regulations, (i.e., Earthing Impedance, Leakage Current).

Instructions for replacing the Controller clock battery

Clock Battery Specifications	
Cell #	CR2032
Classification	Lithium Coin
Chemical System	Lithium / Manganese Dioxide (Li/MnO ₂)
Designation	ANSI / NEDA-5004LC, IEC-CR2032
Nominal Voltage	3.0 Volts
Typical Capacity	220 mAh

1. Turn Controller over.

Diagram 1





Diagram 2



2. To remove the bottom, unscrew pole clamp screw; unscrew six (6) screws shown below. Next, lift the cover.

Diagram 3





Diagram 4



3. Remove power supply assembly
 - a. On the power supply assembly, unscrew two (2) screws shown below (see Diagram 5). Also remove the screws holding the zip ties to the post for the ferrite bead securement.
 - b. Gently lift assembly by holding pole clamp, and prop up with a small block or box (see Diagram 6).

Diagram 5






Diagram 6



4. Remove display panel
 - a. On the display panel, remove the four (4) screws at the end of the panel next to the display screen (see Diagram 7).
 - b. Pull the display panel away from the fan assembly.
 - c. With pliers, gently pull out battery.
 - d. Insert new battery in the same direction as the old battery. Push firmly into place.
 - e. Reinsert and tighten all screws on the display board.




5. Reassembly

- a. Reseat fan assembly and display board (see Diagram 5).
- b. When reseating display board, make sure that all wires are between the two white posts (see Diagram 6).
- c. When reseating the fan assembly on the display board end, make sure that the wires are **inside** the corner post_1 (See Diagram 5) and wrap around the edge of the fan assembly.
- d. On the opposite end from the display board, make sure the wires go **around** the outside of the corner post_2.
- e. Reinsert and tighten the two (2) screws for the power assembly (see Diagram 5).
- f. Reinstall the two screws for the zip ties that were taken apart in step 3a.
- g. Attach bottom cover; reinsert and tighten the six (6) related cover screws (see Diagram 3).
- h. Reinsert pole clamp screw.
- i. Dispose of the old battery in accordance with any government regulations in effect in your area.

Instructions for changing the Controller fuse

Note ⚡ : Remove power cord before inserting screwdriver.

1. Insert a screwdriver on either side of the fuse box and push gently to pop the fuse case out of its socket.



2. Pull the fuse out of the case.
3. Replace the fuse with (Bussmann #S500-5-R or equivalent). The socket in the Controller is keyed so that the fuse case can only be inserted in the correct orientation. Push the case in gently, but firmly, until it snaps into place. At this point, it will be flush with the surrounding surface.

enFlow fluid warming system temperature control and alarms

Temperature control

The enFlow IV fluid/blood warming system includes multiple safety features to prevent overheating of infusion fluids. There are built-in redundancies and back-up safety systems as means for protection in the event of a failure in the primary control mechanism. Several aspects of the system work together to accomplish desired and appropriate safety:

1. Closed-loop temperature control software
2. Audible and visual alarms
3. Software system monitor
4. Independent heater temperature monitoring circuit

Audible/Visual alarms

The enFlow system incorporates an audible/visual alarm system. This system monitors the fluid temperature for an over-temperature condition and system faults. The alarm function is tested each time a Cartridge is inserted. The audible alarm is briefly sounded, and the red light emitting diode indicator is also briefly illuminated. If a dangerous condition occurs, a continuous audible and visual alarm is activated. Over-temperature conditions are calculated according to ASTM Standard F2172-02 (an FDA recognized standard for fluid warmers). This standard allows for spikes in the fluid temperature without activating an alarm. The alarm has a linear time-temperature relationship; therefore, the hotter the fluid the less time it will take for the alarm to be activated. At 45 °C the alarm will sound after about 20 seconds while at 50 °C it will be essentially instantaneous. The alarm will be activated by either the Warmer over-heating the fluid or if the fluid entering the Warmer is too hot. The audible aspect of the alarm can be muted for 1 minute by pressing any key on the Controller. The alarm can also be ended completely by sliding open the covers on the Warmer. (For further information, please reference test for Over-temperature Alarm located in Appendix E: Preventive Maintenance Procedure.)

enFlow troubleshooting

Electromagnetic interference

- ECG, EEG or EMG (cardiac or neuro monitoring) artifact or other interference caused by the enFlow is an uncommon event.
- Cardiac or neuro monitoring interference is common and well-documented in medical literature.
- There are published suggestions to reduce or eliminate the interference that should be employed.
- Also refer to Appendix A Guidance and Manufacturer Declaration - Emissions.

“Interference of the monitored or recorded electrocardiogram is common within operating room and intensive care unit environments.”¹ The enFlow IV fluid/blood warming system, as with all electrical devices, can be associated with some electromagnetic interference (EMI); however, it has been uncommon and inconsistently experienced. Below are troubleshooting suggestions for situations where interference is observed:

Interference confirmation

Turn the unit off. Turn the power supply on the back of the Controller to the OFF position. Reassess the interference. Knowing that cardiac or neuro monitoring is being affected, determine if the interference adversely affects your ability to care for the patient. Consider attempting to reduce the level of interference by employing some simple and readily available solutions.

Interference reduction steps

Check the monitoring pads

All monitoring pads should be fully adhered to the patient's properly prepared skin. Confirm that the pad's foam insulator is not curled up, peeled back or otherwise exposing the conductive gel layer. Confirm that the leadwires' connectors are properly and fully attached. Consider reapplying monitoring electrodes if there is any suspicion they have dried out.

Confirm the patient is properly grounded

In many cases, and in all cases using mono-polar or bi-phase cautery, a grounding pad should be present and applied according to the manufacturer's instructions. As previously suggested, confirm the ground pad is fully adhered to a properly prepared skin surface.

Confirm the enFlow and the ECG monitor are plugged into different outlets

There are two reasons for this action. It is possible that the two systems are in an electrical phase related conflict, which is being expressed on the monitor. Secondly, it is possible that the outlets are not properly grounded or grounded in different locations. While rare, it can be the case and outwardly there would be no way to tell.

Confirm the Warmer cord is not entwined or near the ECG lead cable

Separating the two cords will allow each cable's shielding to work to its full potential.

Confirm the monitoring cables and lead wires are in proper working order

The insulating layer on lead wires and cables degrades over time and with use. Confirm the insulation is intact and operates at its stated specifications.

¹ Patel, Santosh I., M.D., F.R.C.A. and Souter, Michael J, M.B., Ch.B., F.R.C.A.; ***Equipment-related Electrocardiographic Artifacts, Causes, Characteristics, Consequences, and Correction***; Anesthesiology 2008; 108:138–48.

Review the monitor’s notch filter

Check to ensure that the monitoring system’s frequency filter is set appropriately.


Determine the monitor’s sensitivity setting

Many physiological monitors have the ability to interpret electrical signals in two distinct modes: a highly sensitive “Diagnostic” mode or a more filtered “Monitoring” mode. Determine the current mode of operation. If the current mode is set on “Diagnostic” consider adjusting it to “Monitoring.”

Check the ECG pads impedance

Contact your current supplier of monitoring pads or your local Vyair representative to determine if a lower impedance version is available. High impedance monitoring pads are less sensitive to the very low signal strength from the heart beat and appear to be more prone to pronounced interference.

enCheck (Model 400) user guide



Intended use

The enFlow IV Fluid Warmer is indicated for warming blood, blood products and intravenous solutions prior to administration. It is intended to be used by healthcare professionals in hospital, clinical and field environments to help prevent hypothermia.

Purpose

The enCheck tester was developed to quickly and reliably trigger the over-temperature alarm condition on the enFlow Warmer. Within seconds, the enCheck unit will heat the Warmer to an over- temperature scenario causing the alarm to sound. Additionally, the enCheck is designed to verify the Warmer operation at enFlow’s installation site.

Procedure

When the enCheck is connected and running in the normal mode, the heat is generated from the Warmer unit using the same technology as when a cartridge is installed. This mode allows for confirmation of the temperature output of the Warmer. (See Appendix E, section on “Simulated Use Performance Testing.”)

enCheck system description

The enCheck is designed to verify the Warmer operation at enFlow’s installation site. In the Normal Mode, it allows for confirmation of the temperature output of the Warmer. In the Overheat Mode, it heats the Warmer to an over-temperature scenario causing the alarm to sound.

There are three LEDs on the enCheck.

- 1. Lighted Orange LED indicates the power is on.
- 2. Lighted Green LED indicates the enCheck is in the Normal Mode.
- 3. Lighted Red LED indicates the enCheck is in the Overheat mode.

There is a slide switch for switching the overheating on/off.

There is a thermocouple adapter for connecting to a calibrated thermometer for verification.

Environmental Requirements	Performance
Operating Temperature Range: -5 °C to 50 °C	Input Voltage: 12-30 V--- 0.5A
Storage Temperature Range: -30 °C to 70 °C	Maximum of 14 W
Relative Humidity Range: 10% to 90%	

Cleaning the enCheck

1. Use only approved cleaning solutions. (Please reference "Cleaning the enFlow IV fluid/blood warming system components" for a list of approved cleaning solutions.)
2. Moisten a clean cloth with the cleaning solution; do not spray or pour cleaning solutions directly onto the enCheck.
3. Wipe the surface of the enCheck, taking care not to leave excess residual cleaner on the enCheck. If fluid ingress is detected, set the enCheck aside for an extended period of time to allow it to dry.

Appendix A: Technical specifications

Size	Warmer: 12.7 cm L x 6.6 cm W x 4.4 cm H (5.0 in L x 2.6 in W x 1.75 in H) Controller: 23.6 cm L x 16.8 cm W x 9.7 cm H (9.3 in L x 6.6 in W x 3.8 in H) Disposable Cartridge: 11.4 cm L x 3.8 cm W x 1.0 cm H (4.5 in L x 1.5 in W x 0.4 in H) Extension set: 119 mm L x 10.6 mm W (4.7 in L x 0.4 in W)
Weight	Warmer (w/o disposable): 330 g (11.6 oz) Controller: 1.9 kg (4.2 lb) Disposable Cartridge: 33 g (1.2 oz) Extension set: 2 g (0.07 oz)
Parylene Coated Disposable Cartridge (and optional IV extension set):	
Priming Volume	Disposable Cartridge: 4 mL (Optional IV extension set): 0.5 mL
Sterility	Gamma Sterilized
Biocompatibility	ISO 10993
Infusion Set Compatible	ISO 8536-4
Performance:	
Fluid Temperature Output	40 °C ± 2 °C
Flow Rate Range	KVO (2 mL/min) - 200 mL/min
Input Voltage	Warmer: 28.5 VDC at a maximum of 350 Watts Controller: 100-240 VAC
Temperature Set Point	40 °C
Over Temperature Set Point	ASTM F-2172-02
Alarms	IEC 60601-1-8:2006 Ed 2 + A1:2012
Max Input Current	5 A
Input Frequency Range	Warmer: DC Controller: 47-63 Hz
Environmental/Physical Requirements:	
Temperature, Operating	-5 °C to 50 °C
Temperature, Storage	-30 °C to 70 °C
Water Resistance	Warmer: IEC 529 IP67 30 minutes immersion at a depth of 91.4 cm (36 in) Controller: IEC 529 IP31 dripping water Disposable Cartridge (and optional IV extension set): IEC 529 IP68, (IV extension set): IEC 529 IP68 continuous immersion
Penetration	Warmer: IEC 529 IP67 dust tight Controller: IEC 529 IP31 ≥ 2.5 mm diameter Disposable Cartridge (and optional IV extension set): IEC 529 IP68 dust, (IV extension set): IEC 529 IP68 dust tight
Electrical Safety	CAN/CSA C22.2#60601-1:2014 Ed 3, AAMI ES60601-1:2005 + A1:2012, IEC 60601-1:2005 Ed 3 + A1:2012, IEC 60601-1-6:2010 Ed 3 + A1:2013, IEC 60601-1-8:2006 Ed 2 + A1:2012
Software in Medical Devices	IEC 62304:2006 Ed 1 + A1:2015
Relative Humidity, Operating and Storage	Warmer: 10% to 90% Controller: 10% to 90% Disposable Cartridge (and optional IV extension set): 10% to 90%
Altitude, Operating and Storage	Up to 15,000 ft
Air Pressure, Operating and Storage	570 hPa, (17 inHg) to 1060 hPa (31 inHg)
Radiated Magnetic Field Emissions	EMC IEC 60601-1-2:2014
Safety Classifications:	
Type of protection against electrical shock	Class I
Degree of protection against electric shock	Type BF, Defibrillation-Proof
Mode of operation	Continuous

Note : Electromagnetic Compatibility (EMC)

The enFlow IV fluid/blood warming system has been tested and found to comply with the limits for medical devices as set forth in IEC 60601-1-2: (2014) and related standards. These limits are designed to provide reasonable protection against electromagnetic interference (EMI) in a typical medical installation. The enFlow system generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. Harmful interference to other devices can be determined by turning this equipment ON and OFF. Try

to correct the interference using one or more of the following:

* Reorient or relocate the receiving device.

* Increase the separation between the equipment.

* Connect the equipment to an outlet on a different circuit from that to which the other device(s) are connected and consult the Vyair Technician for help.

Care must be taken when operating this equipment around other equipment to avoid reciprocal interference.

Potential electromagnetic interference (EMI), electrostatic discharge (ESD) or other interference could occur to this or to the other equipment and it could interrupt the operation of the device or degrade the performance, which may result in patient injury. Try to minimize this interference by not using other equipment in conjunction with the device.

Guidance and Manufacturer's Declaration – Emissions				
The enFlow is intended for use in the electromagnetic environment specified below. The customer or user of the enFlow should ensure that it is used in such an environment.				
Emissions Test	Compliance		Electromagnetic Environment – Guidance	
RF Emissions CISPR 11	Group 1		The enFlow uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF Emissions CISPR 11	Class A		The enFlow is suitable for use in all establishments, other than domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonics IEC 61000-3-2	Class A			
Flicker IEC 61000-3-3	Complies			
Guidance and Manufacturer's Declaration – Immunity				
The enFlow is intended for use in the electromagnetic environment specified below. The customer or user of the enFlow should ensure that it is used in such an environment.				
Immunity Test	IEC 60601 Test Level	4 th Edition Test Levels	Compliance Level	Electromagnetic Environment- Guidance
ESD IEC 61000-4-2	±6kV Contact ±8kV Air	±8kV Contact ±15kV Air	±8kV Contact ±15kV Air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the r/h should be at least 30%.
EFT IEC 61000-4-4	±2kV Mains ±1kV I/Os	±2kV Mains ±1kV I/Os 100 kHz Repetition Freq	±2kV Mains ±1kV I/Os 100 kHz Repetition Freq	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV Differential ±2kV Common	±1kV Differential ±2kV Common	±1kV Differential ±2kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips/Dropout IEC 61000-4-11	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	0% Ur for 0.5 cycle @0, 45, 90, 135, 180, 225, 270 and 315 degrees. 0% Ur for 1 cycle.	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds ----- 0% Ur for 0.5 cycle @0, 45, 90, 135, 180, 225, 270 and 315 degrees. 0% Ur for 1 cycle.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the enFlow requires continued operation during power mains interruptions, it is recommended that the enFlow be powered from an uninterruptible power supply or battery.
Power Frequency 50/60Hz Magnetic Field IEC 61000-4-8	3 A/m	30A/M 50 OR 60Hz	30 A/M 50/60Hz	Power frequency magnetic fields should be that of a typical commercial or hospital environment.
Guidance and Manufacturer's Declaration – Immunity				
The enFlow is intended for use in the electromagnetic environment specified below. The customer or user of the enFlow should ensure that it is used in such an environment.				

Immunity Test	IEC 60601 Test Level	4 th Edition Test Levels	Compliance Level	Electromagnetic Environment – Guidance		
<div>Conducted RF IEC 61000-4-6</div> <div>Radiated RF IEC 61000-4-3</div>	<div>3 Vrms (Outside ISM) 10 Vrms (In ISM Bands) 150 kHz to 80 MHz</div> <div>10 V/m 80 MHz to 2.5 GHz</div>	<div>3 Vrms (Outside ISM) 6Vrms (In ISM Bands) 150 kHz to 80 MHz</div> <div>10 V/m 80 MHz to 2.7 GHz</div>	<div>(V1) = 3 Vrms (V2) = 6Vrms (V3) = 10Vrms</div> <div>(E1) = 10V/m</div>	<div>Portable and mobile communications equipment should be separated from the enFlow by no less than the distances calculated/listed below:</div> <div>D=(3.5/V1)(Sqrt P) D=(12/V2) (Sqrt P)</div> <div>D=(12/E1)(Sqrt P) 80 to 800 MHz</div> <div>D=(23/E1)(Sqrt P) 800MHz to 2.7 GHz</div> <div>where P is the max power in watts and D is the recommended separation distance in meters.</div> <div>Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1).</div> <div>Interference may occur in the vicinity of equipment containing a transmitter.</div>		
Enclosure Port Immunity to RF Wireless Communications Equipment – IEC 61000-4-3						
Test Frequency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	Immunity Test Level (V/m)
385	380-390	TETRA 400	PM 18Hz	1.8	0.3	27
450	430-470	GMRS 460 FRS 460	FM ±5kHz Dev 1kHz Sine	2	0.3	28
710	704-787	LTE Band 13, 17	PM 217Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5	PM 18Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1, 3, 4, 25 UMITS	PM 217Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE Band 7	PM 217Hz	2	0.3	28
5240	5100-5800	WLAN 802.11a/n	PM 217Hz	0.2	0.3	9
5500						
5785						

Recommended Separation Distances for the enFlow

The enFlow is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the enFlow can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the enFlow as recommended below, according to the maximum output power of the communications equipment.

Max Output Power (Watts)	Separation (m) 150kHz to 80MHz Non-ISM $D=(3.5/V1)(\text{Sqrt } P)$	Separation (m) 150kHz to 80MHz ISM $D=(12/V2)(\text{Sqrt } P)$	Separation (m) 80 to 800MHz $D=(12/E1)(\text{Sqrt } P)$	Separation (m) 800MHz to 2.7GHz $D=(23/E1)(\text{Sqrt } P)$
0.01	.116667	0.12	0.12	0.23
0.1	.368932	0.379473	0.379473	.727324
1	1.166667	1.2	1.2	2.3
10	3.689324	3.794733	3.794733	7.273239
100	11.66667	12	12	23

Sound Levels (with no alarms active)

Sound pressure level Not to exceed 60 dBA at one meter per 60601-1-8

Alarm Sound Level

Alarm	Alarm Priority	Tolerance.
Alarm Volume (Measured per 60601-1-8 Section 6.3.3.2c)	High Priority Alarm: 60dBA	+/- 5dBA

Appendix B: Fault code table

Fault Code Displayed on Controller

#	Name	Description
10	Fluid out temperature sensor reading failure	Sensor is reading a temperature outside the normal operating range or has been determined to be open or shorted.
11	Fluid in temperature sensor reading failure	Sensor is reading a temperature outside the normal operating range or has been determined to be open or shorted.
12	Heater out temperature sensor reading failure	Sensor is reading a temperature outside the normal operating range or has been determined to be open or shorted.
13	Heater in temperature sensor reading failure	Sensor is reading a temperature outside the normal operating range or has been determined to be open or shorted.
20	Average VCC failure	Warmer internal operating voltage is out of range.
21	Average TRef failure	Warmer internal reference voltage is out of range.
30	RAM memory test failure	Data is written to and then read from RAM and verified at power up.
31	ROM memory test failure	A check sum is done on the memory at power up.
32	Flash memory test failure	A check sum is done on the memory at power up.
40	OT fuse or low battery failure	This fault code is indicative of a non-recoverable heater over-temperature, fuse open, condition.
41	Covers failure	This fault occurs every time the covers are opened while heating.
50	Heater MOSFET failure or Circuit Breaker Trip	This fault code is indicative of a Heater Over-temperature protection fault.
52	Clock oscillator failure	This fault is activated if the external crystal oscillator fails.
53	Internal temperature °C failure	This failure is indicative of the Warmer being too hot inside and is caused by operating in an environment that is too warm.
54	Error clock failure	A system clock failure has occurred.
99	General system fault	The Controller is unable to determine the fault.
ACC	Accelerometer fault	Controller accelerometer fault at power up.

Please follow these instructions if any of the fault conditions listed above are displayed on the controller:

1. Thoroughly clean the warmer using the instructions provided in this manual (refer to section titled "Cleaning the enFlow IV fluid/blood warming system components.")
2. Evaluate the exterior of the device. If there is any physical damage to the housings or power cord that would allow fluid ingress into the device, or if the heating element surface is physically damaged, deformed, or breached, then remove the device from service and return the device using the RGA instructions provided in this manual (refer to "Servicing the enFlow IV fluid/blood warming system components.")
3. If there is no obvious physical damage to the device, for Fault Code 50 cycle the power and put back into service. For all other Fault Codes, evaluate the device using the instructions in Appendix E and Appendix F.
4. If the device passes the evaluation, then it may be returned to normal use.

Appendix C: Warming system response by temperature

Fluid Temp	Heater	Temperature LED on Warmer	Display on Controller	Audible Alarm
30 °C	Active	Blue Flashing	30 °C Blue Flashing	No
31 °C	Active	Blue Flashing	31 °C Blue Flashing	No
32 °C	Active	Blue Flashing	32 °C Blue Flashing	No
33 °C	Active	Blue	33 °C Blue	No
34 °C	Active	Blue	34 °C Blue	No
35 °C	Active	Green	35 °C Green	No
36 °C	Active	Green	36 °C Green	No
37 °C	Active	Green	37 °C Green	No
38 °C	Active	Green	38 °C Green	No
39 °C	Active	Green	39 °C Green	No
40 °C	Active	Green	40 °C Green	No
41 °C	Off	Green	41 °C Green	No
42 °C	Off	Green	42 °C Green	No
43 °C	Off	Yellow	43 °C Yellow	No
44 °C	Off	Yellow	44 °C Yellow	No
45 °C	Off	Red Flashing High Priority Alarm after 20 seconds	45 °C after 20 seconds Red Flashing "Over Temp" message	After 20 seconds
46 °C	Off	Red Flashing High Priority Alarm after 16 seconds	46 °C after 16 seconds Red Flashing "Over Temp" message	After 16 seconds
47 °C	Off	Red Flashing High Priority Alarm after 12 seconds	47 °C after 12 seconds Red Flashing "Over Temp" message	After 12 seconds
48 °C	Off	Red Flashing High Priority Alarm after 8 seconds	48 °C after 8 seconds Red Flashing "Over Temp" message	After 8 seconds
49 °C	Off	Red Flashing High Priority Alarm after 4 seconds	49 °C after 4 seconds Red Flashing "Over Temp" message	After 4 seconds
50 °C	Off	Red Flashing High Priority Alarm (immediately)	Red Flashing Over Temp message (immediately)	Immediately

Appendix D: Parts list

Part number	Part	Part number	Instructions For Use
980105VS	Warmer	44000024	System Manual (USA)
980121EU	Controller	44000073	System Manual German (DE)
980200EU	Disposable Cartridge	44000074	System Manual Danish (DA)
980202EU	Disposable Cartridge with IV extension	44000075	System Manual Spanish (ES)
980305VS	Warmer holder	44000076	System Manual Finnish (FI)
980309VS-20	Warmer cord clip	44000077	System Manual French (FR)
980400	enCheck alarm testing tool	44000078	System Manual Italian (IT)
91000178	Power Cord USA	44000079	System Manual Dutch (NL)
91000170	Power Cord Continental Europe	44000080	System Manual Norwegian (NO)
91000172	Power Cord Great Britain	44000081	System Manual Swedish (SV)
91000173	Power Cord Italy	44000083	System Manual Chinese (ZH-S)
91000174	Power Cord Israel	44000084	System Manual Turkish (TR)
91000171	Power Cord Switzerland	44000085	System Manual Portuguese (PT)
91000175	Power Cord India	44000113	System Manual Bulgarian (BG)
91000176	Power Cord Denmark	44000114	System Manual Croatian (HR)
91000177	Power Cord South Africa	44000115	System Manual Czech (CS)
91000179	Power Cord China	44000116	System Manual Greek (EL)
91000180	Power Cord Australia	44000117	System Manual Polish (PL)
91000181	Power Cord New Zealand	44000118	System Manual Russian (RU)
980304VS30	Insulated Strap	44000119	System Manual Serbian (SR)
980331VS-200	Warmer IV Clip	44000123	System Manual Latvian (LV)
980330VS-1	Controller Pole Clamp	44000124	System Manual Lithuanian (LT)
91000178USL	Locking Power Cord- US and Canada	44000125	System Manual Romanian (RO)
91000183	Power Cord Brazil	44000126	System Manual Slovakian (SK)
91000182	Power Cord Japan	44000127	System Manual Slovenian (SL)
		44000128	System Manual Hungarian (HU)
		44000140	System Manual Estonian (ET)
		44000141	System Manual Korean (KO)
		44000142	System Manual Japanese (JA)
		44000143	System Manual Traditional Chinese (ZH-T)

To order the parts below, call the respective manufacturer.

Part number	Part /manufacturer
S500-5-R	Fuse – Bussmann®
CR2032	Battery – Panasonic®

Appendix E: Preventive maintenance procedure

Employ local regulations to determine the frequency of required testing (i.e. Earthing Impedance, Leakage Current) for the enFlow Warmer (Model 100) and Controller (Model 121). Do not repair or perform maintenance while in use with the patient. Before performing maintenance, disconnect the device from the patient, shut it down, and disconnect it from AC Power. All operator maintenance must be performed with the patient off the device. Failure to do so can result in electric shock to the patient and operator.

	Frequency
Functional and Operational Testing Protocols	As required by accrediting body or once a year
Inspections	X
Temperature Readout Display and Status Indicator Lights	X
Electrical Safety	X
Simulated Use Performance Testing: enCheck Model 400 or Alternative Method	X
Alarm Test: enCheck Model 400 or Alternative Method	X

Inspections

1. Ensure all cords and connectors are in good condition and void of any cuts, cracks, or frays. Discoloration from cleaning solutions and disinfectants is normal and to be expected.
2. Ensure the unit is clean and void of any cracks or other signs of damage. If signs of damage are visible, remove it from service and contact Vyaire as soon as possible.

Temperature readout display and status indicator lights

1. Plug the Controller into a functioning power supply. Set the MAINS power to ON. Confirm the Controller power indicator is illuminated and displaying a green color. Confirm the display panel (Controller only) shows in yellow the conditional message "not heating."
2. Connect the Warmer without a Disposable Cartridge inside to the Controller. Confirm the beep signaling connection. Confirm the display continues to show the conditional message "not heating." Confirm the Warmer power LED is flashing green.

Electrical safety

1. Ground wire resistance

Equipment

enFlow Controller
Safety analyzer with test lead

Purpose

The purpose of this test is to check the resistance in ohms of the ground pin to the chassis. For purposes of this check, the pole screw will be considered to be the ground pin, and the chassis is the Controller.

Procedure

A. USA (Tests the equipment inclusive of its power cord.)

1. Remove the Controller from the IV pole.
2. Reinsert the pole clamp screw; screw in tight against the case. (Do not over-tighten.)
3. Attach the banana end of the ground lead on the safety analyzer to the pole clamp screw.
4. Plug the power cord of the Controller into the safety analyzer.
5. Set the function knob on the safety analyzer to ground wire resistance.
6. Set the ground switch to normal.
7. Set the polarity switch to the off position.
8. Power the safety analyzer by plugging in and setting the MAINS power to ON.
9. Record the resistance reading. An acceptable reading is a maximum of 500.0 mΩ.

B. TUV (Tests the equipment exclusive of its power cord.)

1. Follow steps 1-3 in Procedure A.
2. Plug the cord of the safety analyzer into the Controller.
3. Follow steps 5-8 in Procedure A.
4. Record the resistance reading. An acceptable reading is a maximum of 100.0 mΩ at a current of 25 A.

2. Leakage current at the AC power cord

Equipment

Safety analyzer with test lead
enFlow Controller

Purpose

This test is run to check the chassis leakage in microamps. For the purpose of this test, the Controller is the chassis.

Procedure


- 1. Plug the power cord of the Controller into the safety analyzer.
- 2. Turn the knob on the safety analyzer to the chassis leakage function.
- 3. Power the safety analyzer by plugging in and setting the MAINS power to ON.
- 4. Record the polarity and ground readings for the Controller for both power on and off scenarios for the following configurations:

Allowable values of continuous LEAKAGE CURRENTS, in μ A.		
ENCLOSURE LEAKAGE CURRENT	TYPE BF	
Normal polarity — normal ground	Normal Condition	100
Reverse polarity — normal ground	Single Fault Condition	500
Reverse polarity — open ground	Double Fault Condition	500
Normal polarity — open ground	Single Fault Condition	500

3. Leakage current of the Warmer to the saline in the IV line


Equipment

enFlow Warmer and Disposable Cartridge
Safety analyzer with ECG leads
Saline
Wire
Extension set with a non-venting cap




Purpose

The purpose of this test is to check the leakage of current from the Warmer into saline. As IV fluids are generally conductive, a fluid warmer is considered to be electrically connected to the patient similar to an ECG lead; therefore, the leakage needs to be tested. The setup for this procedure is described below.



Procedure

1. Leave the cap on the female end of the Disposable Cartridge.
2. Put an extension set on the male end.
3. Insert a piece of wire into the extension set tubing.
4. Fill the Disposable Cartridge and extension set tubing with saline.
5. Next, put a non-venting cap on the open end of the extension set.
6. Place the Disposable Cartridge setup in the Warmer.
7. Connect the ECG lead from the leakage tester to the wire inserted into the extension set tubing.
8. Perform the ECG lead leakage test.

Note : It doesn't matter which way the wire is inserted into the tubing, just be certain that the wire is in contact with the fluid.


Current	Normal	Single Fault
Earth leakage	5 mA	10 mA
Touch/chassis leakage	100 μ A	500 μ A
Patient leakage	100 μ A	500 μ A
From Table I. Leakage current limits (from IEC 60601-1)		

Simulated use performance testing

A. enCheck Model 400

Equipment

enCheck (includes K type probe)
 Thermal thermocouple meter with $\pm 5^\circ\text{C}$ accuracy
 enFlow Controller (Model 121)
 enFlow Warmer (Model 100)




Purpose

The enCheck Tester was developed to quickly and reliably trigger the over-temperature alarm condition on the enFlow Warmer. Within seconds, the enCheck unit will heat the Warmer to an over- temperature scenario causing the alarm to sound. Additionally, the enCheck is designed to verify the Warmer operation at enFlow's installation site.


Procedure

When the enCheck is connected and running in the normal mode, the heat is generated from the Warmer unit using the same technology as when a cartridge is installed. This mode allows for confirmation of the temperature output of the Warmer.


1. Plug the Controller into a hospital grade outlet.




2. Connect the enCheck to the Controller by inserting the male plug end of the enCheck Hubbell connector into the female receptacle on the front face of the Controller. Push it in so that the plug cover is tight against the receptacle.




3. Next, connect the Warmer to the enCheck by inserting the male plug end of the Warmer into the enCheck female receptacle.




4. Insert the temperature probe connector on the enCheck into a thermometer. Set thermometer to "K" type setting.




5. Insert the enCheck unit into the Warmer. The end of the unit is keyed similar to the Cartridge so it will only fit in the correct orientation. Close the covers.



6. Confirm the enCheck is set to the normal mode.




7. Move the MAINS power switch on the back of the Controller to the ON position. Wait for the thermometer to stabilize, ≈ 30 to 60 sec. assuming all equipment is close to 20 °C.



8. The temperature on the thermometer should be 40 °C \pm 2 °C.

Note : If the temperature in step 8 above is not 40 °C \pm 2 °C, take the Warmer unit out of service. Call Customer Service for an RGA. Please reference the RGA section in “Servicing the enFlow IV fluid/blood warming system components” for further information regarding returns.



Equipment

enFlow system

Power source

Infusion pump capable of maintaining up to 200 mL/min

IV line set

Water bath

Source of distilled water or normal fluid - 0.5 L at $20\text{ }^{\circ}\text{C} \pm 2.0\text{ }^{\circ}\text{C}$

2 extension sets - 22.9 cm (9 in)

Thermometer - capable of measuring $10\text{ }^{\circ}\text{C}$ to $60\text{ }^{\circ}\text{C}$ accurate to $\pm 0.1\text{ }^{\circ}\text{C}$, plus 2 type-K isolated probes


Timer

Waste water container

Graduated cylinder, 100 to 500 mL

Purpose

This test is performed to ascertain that the output fluid temperature of the enFlow system, while using a 22.9 cm (9 in) extension set, is $40\text{ }^{\circ}\text{C} \pm 2.0\text{ }^{\circ}\text{C}$ when the input fluid is $20\text{ }^{\circ}\text{C}$ through the flow rates of 25-100 mL/min (Standard Flows). Additionally, it is run to determine that the rise in fluid temperature is $>16.5\text{ }^{\circ}\text{C}$ when the input fluid is $20\text{ }^{\circ}\text{C}$ utilizing flow rates of 100-200 mL/min (High Flows).



Procedure

Measure input and output temperature of fluid: standard flow

1. Set up the enFlow system for normal operation.
2. Attach an IV line set to a 0.5 liter source of fluid at a temperature of $\approx 20\text{ }^{\circ}\text{C} \pm 2\text{ }^{\circ}\text{C}$. Then run the IV line through a pump capable of maintaining up to 200 mL/min or determine flow rate using the graduated cylinder and timer.
3. Next, attach the IV line to the enFlow system.
4. The temperatures for this test should be measured within 22.9 cm (9 in) of both the input and output connections of the Disposable Cartridge. This

step is done by inserting T connectors in the direct fluid paths, which will accommodate a temperature probe. Connect the temperature probes to a thermometer capable of measuring between 10 °C and 60 °C with ± 0.1 °C accuracy.

5. Prime the IV line setup according to standard IV protocols.
6. Confirm the output end of the extension set empties into the waste water container.
7. Power on the enFlow system, and establish a fluid flow of 100 ± 20 mL/min. Then allow at least 20 seconds for the power-on self-test to complete, the temperature display to read a stable temperature, and the temperature probes to stabilize.
8. Record the input fluid temperature. The acceptable temperature range is $20 \text{ °C} \pm 2.0 \text{ °C}$.
9. Record the output fluid temperature. The acceptable temperature range is $40 \text{ °C} \pm 2.0 \text{ °C}$.
10. It is recommended to repeat steps 1-9 for the flow rate of 60 mL a minute.

Measure input and output temperatures of fluid: high flows

11. Repeat steps 1-7 for the high flow rates of 125 ± 20 mL/min, 175 ± 20 mL/min, and 200 mL/min. However, in place of steps 8 and 9, measure the **rise** in temperature of the output fluid over the input fluid value. The rise should be $>16.5 \text{ °C}$.

Alarm test

A. enCheck Model 400

1. enCheck Model 400 (includes "K" type probe)
2. enFlow Controller/AC Power Pack (Model 121)
3. enFlow Warmer (Model 100)
4. Calibrated Thermometer with $\pm 0.5 \text{ °C}$ accuracy (Fluke 54 or equivalent)

Test procedure

Alarm function test:

1. Plug the enFlow Controller into an AC outlet.
2. Confirm that the switch on the enCheck is set to the normal mode.
3. Connect the enCheck to the enFlow Controller, by inserting the male plug end of the enCheck Hubbell connector into the female receptacle on the front face of the enFlow Controller. Push in and twist the enCheck Hubbell connector so that the plug cover is locked snugly against the receptacle.
4. Connect the enFlow Warmer to the enCheck by inserting the male plug end of the Warmer into the female receptacle of the enCheck. Push the plugs together and twist to lock them snugly against each other.
5. Insert the enCheck temperature probe connector into the Calibrated Thermometer and set the Thermometer to "K" type setting.
6. Insert the enCheck unit into the Warmer so the bottom end of the unit is keyed similarly as the Warmer Cartridge, so that it only fits in the correct orientation (curved section on the left side), and close the Warmer covers.
7. Switch the Main power switch on the back of enFlow Controller to the ON or **I** position.
8. Verify that the Orange "power" LED indicator is illuminated on the enCheck.
9. After waiting, for approximately 30 to 60 seconds, confirm the temperature on the Calibrated Thermometer is reading in the range of $40 \text{ °C} \pm 2 \text{ °C}$, and that the Green "Normal" LED on the enCheck is illuminated.
10. Disconnect the enCheck temperature probe connector that is connected to the Calibrated Thermometer.
11. Place the Switch on the enCheck to the "Overheat" position and verify that the Red "Overheat" LED is illuminated on the enCheck.
12. At this time, observe the enFlow Controller display until the temperature rises to greater than 45 °C ; the controller display should be Red and the audible alarm should sound, the Green "TEMP" LED on the enFlow Warmer should also have changed from a solid Green to a flashing Red.
13. Place the Switch on the enCheck to the "Normal" position and re-connect the enCheck temperature probe connector to the Calibrated Thermometer.
14. Verify the LED and display indicators are red. Monitor the temperature on the Calibrated Temperature meter and on the enFlow Controller display; the temperature will gradually drop, and as it falls below 45 °C , the Controller display and the Warmer LED should change state from Red to Yellow. Verify that the Warmer LED and the Controller display indicators change to yellow as the temperature transitions through the 44 °C and 43 °C temperature window; the Controller indicators and the LED will return to green as soon as the temperature drops to 42 °C and will remain green, until the temperature drops to below 35 °C .

Equipment

enFlow system

Power source

Infusion pump capable of maintaining up to 100 mL/min

IV line set

Water bath

Source of distilled water or normal fluid – 0.5 L at $50 \text{ °C} \pm 2.0 \text{ °C}$

2 extension sets - 22.9 cm (9 in)

Thermometer - capable of measuring 10 °C to 60 °C accurate to $\pm 0.1 \text{ °C}$, plus 2 external probes



Waste water container

Purpose

The purpose of this test is to determine that the over-temperature alarm on the Warmer is working properly.

Procedure

1. Use the Simulated Use Performance Testing setup described in steps 1-5 above.
2. Use a clamp to stop the flow in the IV line.



3. Change the source fluid temperature to $50\text{ °C} \pm 2\text{ °C}$. Take extra care while working with hot fluids.
4. Power on the enFlow system.
5. Allow at least 20 seconds for the power-on self-test to complete and the temperature display to reach a stable temperature. Confirm all LEDs are illuminated green.
6. Release the clamp and establish a fluid flow of $100 \pm 20\text{ mL/min}$.
7. With the thermometer, confirm the temperature of the fluid at both the input and output end of cartridge.
8. The Over-Temperature alarm sounds within approximately 20 seconds or less of the output fluid temperature reaching that of the input fluid level.
 - A. Confirm the Controller produces an audible beep and displays a red Over-Temp message on the display (Controller only).
 - B. Confirm the temperature LED on the Warmer flashes red indicating that the fluid has gone over-temperature.

Preventative Maintenance Checklist

[illegible]


Appendix F: enFlow IV fluid/blood warming system operational checklist

Warmer Serial No _____ Controller Serial No _____

Warming system location/identifier _____

Date _____

Procedure Instructions	Pass/Fail	Input Temp	Output Temp
Inspection			
Ensure all cords and connectors are in good condition and void of any cuts, cracks, or frays.			
Ensure units are clean and void of any cracks or other signs of damage.			
Performance Test Setup			
Set up the system for normal operation. Provide a 0.5 liter source of fluid at 20 °C ± 2 °C. Measure the temperature within 22.9 cm (9 in) of both the input and output connections of the Disposable Cartridge by inserting a T connector in the direct fluid path, which will accommodate a temperature probe. Connect the temperature probes to a meter capable of measuring between 10 °C and 60 °C with 0.1 °C accuracy. Prime the IV line setup according to standard IV protocols. Turn the enFlow system on and establish a fluid flow of 100 ± 20 mL/min. Wait for the temperature probes to stabilize.			
Record the input fluid temperature. Input fluid temperature 20 °C ± 2 °C.			
Record the output fluid temperature. Output fluid temperature 40 °C ± 2 °C.			
Over-Temperature Alarm Check			
Use performance testing setup. Change the source of fluid's temperature to 50 °C ± 2 °C. Turn the enFlow system on and establish a fluid flow of 100 ± 20 mL/min. Wait for the temperature at the probes to stabilize.			
Record the input and output fluid temperatures.			
The High Priority Over Temp Alarm occurs within less than 20 seconds of reaching input temperature. (See Appendix C.)			
High Priority Alarm indicated by Audible beep and Over-Temp message in Red appearing on the Controller			
Red Temperature LED flashes on the Warmer, also indicating a High Priority Alarm			
Electrical Safety			
Follow safety analyzer manufacturer's instructions			
Test leakage current at the AC power cord using a safety analyzer.			
Test leakage current of the Warmer to the saline in the IV line using a safety analyzer.			
Inspected By			
Enter initials and confirm date.			
Comments, Observations or Corrective Actions			

Note : Please reference the Preventative Maintenance Procedure (Appendix E) for the specific procedures in order to perform the tests listed above.


Appendix G: enFlow IV fluid/blood warming system operational checklist—enCheck Model 400 method

Warmer Serial No _____ Controller Serial No _____

Warming system location/identifier _____

Date _____

Procedure Instructions	Pass/Fail	Temperature
Inspection		
Ensure all cords and connectors are in good condition and void of any cuts, cracks, or frays.		
Ensure units are clean and void of any cracks or other signs of damage.		
Performance Test Setup		
Plug the Controller into a hospital grade outlet. Connect the enCheck to the Controller by inserting the male plug end of the enCheck Hubbell connector into the female receptacle on the front face of the Controller. Push it in so that the plug cover is tight against the receptacle. Insert the enCheck unit into the Warmer. The end of the unit is keyed similar to the Cartridge so it will only fit in the correct orientation. Close the covers. The temperature on the thermometer should be $40^{\circ}\text{C} \pm 2^{\circ}\text{C}$.		
Record the temperature.		
Over-Temperature Alarm Check		
Plug the enFlow Controller into an AC outlet. Confirm that the switch on the enCheck is set to the normal mode. Connect the enCheck to the enFlow Controller, by inserting the male plug end of the enCheck Hubbell connector into the female receptacle on the front face of the enFlow Controller. Push in and twist the enCheck Hubbell connector so that the plug cover is locked snugly against the receptacle. Connect the enFlow Warmer to the enCheck by inserting the male plug end of the Warmer into the female receptacle of the enCheck. Push the plugs together and twist to lock them snugly against each other. Insert the enCheck temperature probe connector into the Calibrated Thermometer and set the Thermometer to "K" type setting. Insert the enCheck unit into the Warmer so the bottom end of the unit is keyed similarly as the Warmer Cartridge, so that it only fits in the correct orientation (curved section on the left side), and close the Warmer covers. Switch the Main power switch on the back of enFlow Controller to the ON or I position. Verify that the Orange "power" LED indicator is illuminated on the enCheck. After waiting, for approximately 30 to 60 seconds, confirm the temperature on the Calibrated Thermometer is reading in the range of $40^{\circ}\text{C} \pm 2^{\circ}\text{C}$, and that the Green "Normal" LED on the enCheck is illuminated. Disconnect the enCheck temperature probe connector that is connected to the Calibrated Thermometer. Place the Switch on the enCheck to the "Overheat" position and verify that the Red "Overheat" LED is illuminated on the enCheck. At this time, observe the enFlow Controller display until the temperature rises to greater than 45°C ; the controller display should be Red and the audible alarm should sound, the Green "TEMP" LED on the enFlow Warmer should also have changed from a solid Green to a flashing Red.		
Record the temperature.		
The High Priority Over Temp Alarm occurs within less than 20 seconds of reaching input temperature. (See Appendix C.)		
High Priority Alarm indicated by Audible beep and Over-Temp message in Red appearing on the Controller		
Red Temperature LED flashes on the Warmer, also indicating a High Priority Alarm		
Electrical Safety		
Follow safety analyzer manufacturer's instructions		
Test leakage current at the AC power cord using a safety analyzer.		
Test leakage current of the Warmer to the saline in the IV line using a safety analyzer.		
Inspected By		
Enter initials and confirm date.		
Comments, Observations or Corrective Actions		

Note : Please reference the Preventative Maintenance Procedure (Appendix E) for the specific procedures in order to perform the tests listed above.

Appendix H: Glossary

enFlow IV fluid/blood warming system	The enFlow IV fluid/blood warming system consists of three products: the Warmer (No. 980105VS), the Controller (No. 980121EU), and the Disposable Cartridge (No. 980202EU/980200EU), which together form a system designed to warm intravenous fluids and blood products to reduce hypothermia.
Warmer (Model 100)	The Warmer is a small, lightweight, robust fluid warmer that heats blood, blood products, and intravenous fluids to 40 °C from flow rates of KVO to 200 mL/min.
Controller (Model 121)	The Controller displays a temperature readout in degrees C, and features a keypad, to control the clock and the mute feature. Additionally, it converts AC line power to 28.5 Volts DC, and is used as a power source for the Warmer.
Disposable Cartridge (Model 200)	The sterile, single-patient-use Disposable Cartridge is used as an in-line component of an IV infusion set to heat fluids/blood being administered to the patient's body.
Disposable Cartridge with IV extension set (Model 202)	The Disposable Cartridge with IV extension set contains the same Disposable Cartridge described above. In addition, it includes a sterile, single-patient-use IV extension set.
Intravenous fluids	Fluids such as Normal Saline, Dextrose, Dextron, Packed RBC's
enCheck Model 400	enCheck alarm testing tool
KVO	Keep Vein Open refers to an intravenous infusion rate defined as approximately 2 mL/min (121 mL/hr).
LED	Light Emitting Diode
mL/min	Milliliters per minute
RBCs	Packed Red Blood Cells