

LARYNGOSCOPES – SUNMED PORTFOLIO




INTRO TO LARYNGOSCOPY


Laryngoscopy (la·ryn·go·sco·pi):

Laryngoscopy is a term describing the visualization or examination of the larynx by distraction of the upper airway structures.


It is typically used to enable intubation and airway management in anesthesia, critical care, and trauma scenarios.



AIRWAY ANATOMY





AIRWAY VIEW IN SUPINE POSITION



 alamy stock photo

ADTXGE
www.alamy.com

KEY AIRWAY ANATOMY FOR LARYNGOSCOPY:



Glottic opening: the opening between the vocal cords. This is the target for tracheal intubation.

Vallecule (val·lec·u·la): the space between the base of the tongue and the epiglottis

Arytenoid cartilages (ar·y·te·noid): A pair of small triangular cartilages in the larynx that help to move the vocal cords

INTRODUCTION TO LARYNGOSCOPES

A **laryngoscope** (la·ryn·go·scope) is used to lift the upper airway structures, such as the epiglottis, out of the way to allow visualization of the vocal cords (larynx) and enable intubation through the glottic opening.

Laryngoscopes are used where-ever tracheal intubations are performed, including OR, ED, ICU, NICU, PICU, EMS



OPTIONAL VIDEOS

Good views from practitioner's POV

https://www.youtube.com/watch?v=4V_poulbcnA

Direct Laryngoscopy <https://www.youtube.com/watch?v=AZeBumPaj4g>

Use of Mac and Miller blades for intubation

Anatomy and Intubation


<https://youtu.be/iOPpSGbuYmQ>

- Skip to the 6:30 min mark: @ 6:30 Anatomy; Great view of vallecula lift at 7:30 / Mallampati overview
- @ 8:30 Intubation Equipment; @ 12:00 Direct laryngoscopic intubation
The first 6.5 mins are pharmacology – double back to it if you can

Intubation using Direct and Indirect (Video) Laryngoscopes

<https://www.youtube.com/watch?v=gnkYGRMaw7o>

First 6 min shows use of different direct and video laryngoscopy devices




OBJECTIVES & CONTENT OUTLINE

**You will have the knowledge, tools, support and target list to be successful
with IntuBrite Day 1 post-NSM**

Session	Deliverables
I. Laryngoscopy Intro	Review of Key Basics Market Overview
II. Product Details	IntuBrite Products Competitive Products
III. Sales Process	Clear sales process <ul style="list-style-type: none">• Targeting• Detailed Talk Track• Evaluation planning & management• FAQs/ Objection Handling
IV. Planning Your Success	<ul style="list-style-type: none">• Intubation training• Leading an Effective In-Service

WHY IS INTUBRITE SO EXCITING?



INTUBRITE®


INTUBRITE DEMO KITS

Demo Kit Contents:

- **Dual LED Disposable** Handles & Blades
- **Greenlight Disposable** Handles & Blades
- **VLS Edge Video** Laryngoscope, Reusable Blades, Wands, Rigid Styles, Disposable Sheaths, Accessories
- **Competitive Disposable Samples:** Flexicare, Teleflex Rusch
- **Misc accessories:** Cases, PFHV & STLT, cloth

POP QUIZ

Name the Airway Anatomy



- *Glottic opening*: the opening between the vocal cords. This is the target for tracheal intubation.
- *Vallecula* (val·lec·u·la): the space between the base of the tongue and the epiglottis
- *Arytenoid cartilages* (ar·y·te·noid): A pair of small triangular cartilages in the larynx that help to move the vocal cords

POP QUIZ:

Which Blade is Mac Style & which is Miller Style?

A.



Mac Blade

- The most common **Curved** blade design
- **MaC**

B.



Miller Blade

- Most common **straight** blade design
- **MiLLer**


BLADE SIZES – MAC, MILLER SIZE RANGES

- Mac and Miller blades overlap in many, but not all size offerings
- Straight blades are preferred for NICU use – reflected in size range of Miller blades.


000, 00 and 1.5 are important differentiators to NICU and children's facilities

Size	Patient Size	% of market	Mac	Miller
000	Small premature infant	Minimal		✓
00	Premature infant	1.2%		✓
0	Neonate	2.8%	✓	✓
1	Small Child	3.7%	✓	✓
1.5	Child	1.2%		✓
2	Child	19.2%	✓	✓
3	Adult	37.5%	✓	✓
3.5	Adult	1.8%	✓	
4	Large adult	22.6%	✓	✓
5	Extra-large adult	Minimal	✓	

MORE BLADE SHAPES



BLADE DESIGN IMPACTS TECHNIQUE




Straight blade: Tip goes *under the epiglottis* and lifts it directly

Curved blade: Tip fits *into the vallecula* and indirectly lifts the epiglottis


FIBER OPTIC (OR GREENLINE) LIGHTING SYSTEMS

- >90% of market
- Handles and blades built to the Green ISO Standard 7376 can use the term Green, Green System
- Indicated by green color coding on blade and handle
- Why is this important?
 - Indicates compatibility with any other handle or blade that is also “Green”
 - Can mix & match disposable & reusable



STANDARD LED LIGHTING SYSTEMS

- <10% of market
- Generally have a superior lighting system than fiberoptics
fiber optic pipe degrades light quality
- Not compatible with Green systems



POP QUIZ


Which illustrates use of Direct Laryngoscopy ____?

A.



Direct Laryngoscopy

B.



**Indirect, or Video
Laryngoscopy**


LARYNGOSCOPY MARKET



MAJOR COMPETITOR PRODUCT OFFERINGS

	Direct				Video	
	Standard LED or Proprietary		Greenline (FiberOptic LED)		Handheld	Rollstand
	Reusable	Disposable	Reusable	Disposable		
IntuBrite	✓ Handle	✓		✓	✓	✓
Flexicare		✓ (new)		✓		
Teleflex Rusch	✓	✓	✓	✓	✓	
SunMed	✓	✓ Blade	✓	✓		
Vital Signs				✓		
Verathon	✓ Handle	✓ Blade		✓	✓ (new)	✓
Medtronic					✓	
Ambu					✓	
Karl Storz			✓	✓	✓	✓

SALTER IS POSITIONED TO TAKE ADVANTAGE OF LARYNGOSCOPE MARKET TRENDS



OUR GREATEST OPPORTUNITY IS HOSPITALS WITH DISPOSABLE HANDLES AND BLADES

Why Hospitals?

- Perform 19 out of 20 intubations annually
- Leverages our existing relationships
- Leverages our acute product portfolio

Why Disposable?

- 3x the unit growth rate of reusables
 - Driven by JCAHO guidelines
 - Increased focus on HAI
- Ongoing, consistent purchase, not a one time or infrequent purchase

REUSABLE → DISPOSABLE: JCAHO DRIVES CHANGE

FDA and JCAHO follow the Spaulding Classification...			...which results in citations when hospitals are only doing low-level disinfection on handles	
Type	Contact	Req		
Critical	Enters sterile tissue or blood flow	Sterilization		
Semi-Critical	Contacts mucus membrane	High-level disinfection		
Non-critical	No- or intact-skin contact	Low-level disinfection		

— **Correct Practice:** Typically, submerged in solution by central reprocessing

— **Current Practice:** Wipe down between cases

↑ JCAHO citation



Guideline for Disinfection and Sterilization in Healthcare Facilities, 2008

Semicritical Items

Semicritical items contact mucous membranes or nonintact skin. This category includes respiratory therapy and anesthesia equipment, some endoscopes, laryngoscope blades²⁴, esophageal manometry probes, cystoscopes²⁵, anorectal manometry catheters, and diaphragm fitting rings. These medical devices should be free from all microorganisms; however, small numbers of bacterial spores are permissible. Intact mucous membranes, such as those of the lungs and the gastrointestinal tract, generally are resistant to infection by common bacterial spores but susceptible to other organisms, such

JCAHO GUIDELINES

Laryngoscopes Blades and Handles - How to clean, disinfect and store these device


Laryngoscopes – Blades and Handles - How should we clean, disinfect and store these devices? How will the surveyors evaluate this process?

- **Laryngoscope blades** are processed via either high-level disinfection or sterilization.
- **Laryngoscope handles**, the organization is following the manufacturer's instructions-for-use for cleaning/disinfection guidance.
- **Packaged**
- **Stored in a way to prevent contamination**

https://www.jointcommission.org/standards_information/jcfaqdetails.aspx?StandardsFaqlId=1201&ProgramId=46

STERILIZATION STEPS & COST

- High cost per intubation: \$17-\$26
- Laryngoscope failure due to wear & tear, improper re-assembly
- Risk of cross-contamination resulting from poor cleaning or bacteria that remains
- Loss of product , replacement
- Time delays in product availability
- Diminished light quality over time



From 'Flexicare Cost of Reprocessing Brochure'

<http://www.briteprosolo.com/images/downloads/Cost%20of%20Reprocessing.pdf>

Fiber optic laryngoscope handles Directions for use

2 Directions for use

Welch Allyn rechargeable fiber optic and standard laryngoscope handles

storing the device.

Reprocessing instructions

These reprocessing instructions refer to procedures for cleaning and intermediate level disinfection. Rechargeable laryngoscope handles must be reprocessed prior to first use and between each use using the following method as outlined in this document:

Directions for use

Welch Allyn rechargeable fiber optic and standard laryngoscope handles 3

Initial cleaning and disinfection:

1. Follow the germicidal wipe manufacturer's instructions to clean all exposed surfaces of handle parts, bottom cap (and lamp cartridge for 60713 & 60835).
2. If necessary, brush with a dry, soft-bristled brush and re-wipe to loosen/remove excessive visible soil.

Reprocessing instructions

These *reprocessing* instructions refer to procedures for cleaning and intermediate level disinfection. Rechargeable laryngoscope handles must be reprocessed prior to first use and between each use using the following method as outlined in this document:

- **Cleaning and intermediate level disinfection**

Welch Allyn has validated the above instruction as being capable of preparing these laryngoscope handles for re-use. The user must ensure that the reprocessing as actually performed by the user's personnel, with the user's equipment and materials, achieves the desired result. This may require validation and routine monitoring of the user's actual process

ACCESSED 9/20/2017
[HTTPS://WWW.WELCHALLYN.COM/CONTENT/DAM/WELCHALLYN/DOCUMENTS/SAP-DOCUMENTS/LIT/80021/80021285LITPDF.PDF](https://www.welchallyn.com/content/dam/welchallyn/documents/sap-documents/lit/80021/80021285litpdf.pdf)



Figure 1



Figure 2



Figure 3

disinfection

disinfection as directed by the germicidal wipe manufacturer.

Storage

Store handle per facility practice to allow device to remain clean, dry, and ready for service.

End of reprocessing instructions for cleaning and intermediate level disinfection.



REF 901038 LARYNGOSCOPE
DIR 80021285 Ver. A
© 2016 Welch Allyn, Inc.
Revision date: 2016-02

R_x ONLY

Fiber optic laryngoscope handles Directions for use

ACCESSED 9/20/2017
[HTTPS://WWW.WELCHALLYN.COM/C
ONTENT/DAM/WELCHALLYN/DOCUM
ENTS/SAP-
DOCUMENTS/LIT/80021/80021285LITPD
F.PDF](https://www.welchallyn.com/content/dam/welchallyn/documents/sap-documents/lit/80021/80021285litpdf.pdf)

Clean and Intermediate Level Disinfection:

Prepare for Cleaning

- Separate blade from handle
- Prevent handle from drying

Initial Cleaning

- Remove the battery
- Clean with germicidal wipe
- Brush to remove any visible soil
- Re-wipe with germicidal wipe and let sit per germicidal wipe instructions

Air dry

Inspect parts and reassemble

Test by attaching to a clean and disinfected blade

Re-wipe all surfaces with germicidal wipe and let sit per germicidal wipe instructions

Store per facility practice




Figure 3

SIMPLEST WAY TO KEEP JCAHO COMPLAINT


“We use disposable handles and blades”

DIRECT LARYNGOSCOPY → VIDEO LARYNGOSCOPY

Drivers of
conversion to
Video

1. Clinical advantage of video
 - More compelling for less skilled intubationists
 - Strong for EMS, ED, ICU
 - Less for anesthesiologists
2. Clinical advantage for difficult cases
3. Seen as best, newest technology
4. More affordable technology options

VIDEO LARYNGOSCOPY COMPETITION FOR 'ROUTINE'



MAJOR PLAYERS BY SEGMENT

Reusable	
Vyaire	19%
Teleflex	16%
Sun-Med	10%
Tri-Anim	8%
Welch-Allyn	7%
Heine	6%
Propper	5%
Flexicare	5%
All Others Armstrong, Mercury, Medline, Sharn, Anesthesia Service, > 30 others	24%

Fully Disposable	
Flexicare	65%
Teleflex	21%
IntuBrite	12%
All Others (Heine, Welch Allyn, OPB, Sun-Med, Tri-Anim, Vyaire)	2%

Video	
Verathon	#1
Karl-Storz	#2
Medtronic	#3
Ambu	#4
Others Teleflex, Venner, McGaw, DRE, Pentax	Minimal

* Reusable includes disposable blades used with reusable handles

REUSABILITY


Reusable

Reusable items are designed to withstand reprocessing (high-level disinfection or sterilization). Most are guaranteed for duration (e.g. 5-7 years) and/or a defined number of autoclave cycles (e.g.) 4,000. Often, the electronic parts (bulb, battery, etc.) can be replaced when broken

Disposable

Disposable items are designed for single-patient use.

Intubations and Sales by Product Type



Key Takeaways

- All the money is in disposables.
- Winning disposable blades only to use with current reusable handles is a pot. wedge
- Winning reusables is unlikely. Don't devote time here.

REUSABLE VS DISPOSABLE: TRADE OFFS

	Reusable Handle & Blade	Disposable Handle & Blade
Benefits	<ul style="list-style-type: none">• Reusable Handle & Blade• High quality <u>weight and feel</u>• What most physicians are accustomed to using – viewed as the gold standard	<ul style="list-style-type: none">• Fresh battery and light for each use• No maintenance or reprocessing resources• Single patient use – no risk of cross-contamination
Drawbacks	<ul style="list-style-type: none">• Performance failure/variability due to reprocessing wear & tear• Light quality degrades/ flicker• Battery failure, incorrect reinstallation• Component failure• High initial capital investment• Maintenance & replacement• Reprocessing time & costs (chemicals, time, labor, supplies)	<ul style="list-style-type: none">• Not perceived to have the quality of a reusable• Cost is more ‘tangible’ to department• Environmental impact perceived to be more than reusable

REUSABLE → DISPOSABLE: COMPLIANCE CHALLENGES

Switching to high-level disinfection of handles is expensive

\$

Cost of high-level disinfection

- Cleaning solution and sterile water
- Labor costs

\$

Must invest in many more reusable scopes

- Go from 1 set per OR to multiple sets per OR to cover time while being cleaned
- Higher breakage and loss rates during cleaning

Customers are reluctant to believe this

**\$17-\$26 per
intubation for
reprocessing
handle and
blade**

Compliance to high-level disinfection is impractical

- Staff are frequently non-compliant
- Unanticipated breakage due to
 - Leaving batteries in handle during cleaning or incorrect insertion
 - Losing scopes during transport
 - Faster wear-and-tear due to reprocessing

**Going
disposable is the
simplest way to
stay compliant**

COST OF REPROCESSING REUSABLE LARYNGOSCOPES

Outcome

As a result of the study, it was found that Glendale Adventist Hospital is incurring \$129,854 in reprocessing costs directly associated with reusable laryngoscope blades and handles annually. With an estimated 5,000 intubations per year, this results in an average cost/intubation of \$25.97. Table 1 on the following page shows a breakdown of the costs found within this study.

Furthermore, this does not include the potential impact from reducing the occurrence of hospital acquired infections. In addition to patient care issues, this financial liability has been conservatively estimated at \$40,000 per hospital acquired infection occurrence. Assuming that even one reduction in hospital acquired infection is prevented by utilizing single-use laryngoscopes, this offers an additional savings of \$8 per intubation for Glendale Adventist Hospital.

Glendale Adventist Hospital Glendale, CA

- Reprocessing cost \$25.97/intubation
- If 1 HAI: Add'l \$8/intubation

Total Reusable Costs

Yearly costs of sterilization in CS (labor & supplies)

Yearly costs of cleaning in Anesthesia room and OR Checking(labor & supplies)

Purchase cost of new handles and blades

	Re-Usable
Yearly costs of sterilization in CS (labor & supplies)	\$ 44,035
Yearly costs of cleaning in Anesthesia room and OR Checking(labor & supplies)	\$ 14,826
Purchase cost of new handles and blades	\$ 7,688
Total Annual Cost	\$ 66,548
Cost/Intubation	\$ 17.12

Cost of Reusable Laryngoscope System if incurring only 1 Hospital Acquired Infection Per Year

Estimated Cost of One Hospital Acquired Infection


\$ 40,000

Total Annual Cost	\$ 106,548
Cost/Intubation	\$ 27.40



Southeast Health Cape Girardeau, Missouri

- Reprocessing cost \$17.12/intubation
- If 1 HAI: Add'l \$10.28/intubation

MARKET SEGMENT SIZES AND GROWTH






MARKET SIZE



Salter's greatest opportunity is in the Hospital market with
Disposable handles & blades

* Reusable includes disposable blades used with reusable handles

MAJOR PRODUCTS, PRICING, STRENGTHS AND DRAWBACKS: DISPOSABLE

		ASP	Strengths	Weaknesses
	Flexicare BritePro	\$10-15	<ul style="list-style-type: none">• Recycling program• Contracts	<ul style="list-style-type: none">• Reliability• Handle-blade give
	Teleflex DispoLED	\$10-15	<ul style="list-style-type: none">• None	<ul style="list-style-type: none">• Reliability• Light quality• Handle-blade give
	Intubrite	\$11-12	<ul style="list-style-type: none">• UV light• Light quality• Metal construction• Reliability, strength, and quality	

Smaller Players (not comprehensive)

	Heine XP		Sunmed SunOne		Storz Laryngobloc		OBP Surescope
---	-----------------	---	----------------------	---	--------------------------	---	----------------------

MAJOR PRODUCTS AND PRICING: VIDEO

				Reusable Config				
			Monitor	Blades	Cable	Wands	Blades	
ROLLSTAND		Verathon Glidescope	~\$5k	\$2-3k	~\$5-\$6k	~\$4k	~\$20-30	
		Karl-Storz C-MAC	Over \$5k	~\$4k	~\$500	unknown	unknown	
		Intubrite VLS 8800	\$4-\$5k	\$1k	N/A	~\$1k	\$8-10	
HANDHELD		Ambu King Vision	\$1.0-\$1.5k	N/A	N/A	included	\$8-10	
		Medtronic McGrath MAC	\$1.6-\$2.6k Sometimes free	N/A	N/A	included	\$12-20	
		Intubrite 6600 Edge	\$1.5 - \$2.5k	\$1k		~\$1k	\$8-10	
Smaller Players (not comprehensive)		Venner APA dist. by Vyaire		Airtraq Dist. by TFX		Magaw Co-Pilot	 Karl Storz C-MAC PM	Verathon Glidescope Go (NEW)

QUESTIONS?



 SunMed