

PHYSIOLOGY AND ANATOMY LARYNGOSCOPES




INTRO TO LARYNGOSCOPY


Laryngoscopy (la·ryn·go·sco·pi):

Laryngoscopy is a term describing the visualization or examination of the larynx by distraction of the upper airway structures.


It is typically used to enable intubation and airway management in anesthesia, critical care, and trauma scenarios.



AIRWAY ANATOMY





AIRWAY VIEW IN SUPINE POSITION



 alamy stock photo

ADTXGE
www.alamy.com

KEY AIRWAY ANATOMY FOR LARYNGOSCOPY:



Glottic opening: the opening between the vocal cords. This is the target for tracheal intubation.

Vallecule (val·lec·u·la): the space between the base of the tongue and the epiglottis

Arytenoid cartilages (ar·y·te·noid): A pair of small triangular cartilages in the larynx that help to move the vocal cords

INTRODUCTION TO LARYNGOSCOPES

A **laryngoscope** (la·ryn·go·scope) is used to lift the upper airway structures, such as the epiglottis, out of the way to allow visualization of the vocal cords (larynx) and enable intubation through the glottic opening.

Laryngoscopes are used where-ever tracheal intubations are performed, including OR, ED, ICU, NICU, PICU, EMS



DIRECT VS. INDIRECT (VIDEO) LARYNGOSCOPY


Direct Laryngoscopy:

- Uses direct vision or line of site from the eye of the clinician to the patients tracheal opening or glottic opening.



Indirect Laryngoscopy


- Techniques that do not use direct glottic opening visualization. In video laryngoscopy, an image of the patient glottic opening is transmitted to a monitor screen.



LARYNGOSCOPY USE CASES FOR TRACHEAL INTUBATION

Est. % of Intubations	Dept	Primarily used by	Use Environment	Use profile
75%	OR	Anesthesiologist CRNA	Planned, Controlled	<ul style="list-style-type: none"> • Highest volume of placements • Largely direct laryngoscopy • Limited video laryngoscopes, primarily used for difficult airway
7%	ICU NICU Floor	RT Intensivist	Emergent	<ul style="list-style-type: none"> • Laryngoscopes in crash cart • Use is infrequent & unplanned
12%	ED	Intensivist	Emergent	<ul style="list-style-type: none"> • Often difficult, traumatic airway
5%	EMS	Paramedic EMT	Emergent, uncontrolled environment	<ul style="list-style-type: none"> • Highest use of video • Difficult, traumatic airway • Less frequent intubations so less experienced in laryngoscopy


LARYNGOSCOPES ARE DEFINED BY SEVERAL VARIABLES



BLADES TYPES - INTRODUCTION

There are numerous blade designs, but they have common features:

- Spatula – curved surface that passes over the surface of the tongue
- Flange – used to direct or displace the tongue. The shape and height of flange can vary significantly between blade designs and manufacturers
- Tip – designed to lift the epiglottis either directly or indirectly
- Light source – provides illumination of the anatomy and vocal cords



BLADE SHAPES: CURVED VS. STRAIGHT

- All blades generally fall into either a “Curved” or “Straight” style.

Curved



- Wider flange and curved blades are helpful in keeping the tongue retracted from the field of vision and allowing more room for passing the tube

Straight




- Provide straighter path for the endotracheal tube, making intubation easier once the larynx can be seen.
- Preferred in patients with large, floppy epiglottis because the tip is designed to lift the epiglottis out of the way
- More likely to be used in pediatrics


THE MOST COMMON BLADES: MAC AND MILLER

Although there are numerous blade designs, more than 90% of intubations are done with either a Macintosh, aka Mac, or Miller blade.

Mac Blades
(Curved, like the 'C' in Mac)



Miller Blades
(Straight, like the 'L's' in Miller)



BLADE SIZES


Size	Patient Size	% of market
000	Small premature infant	Minimal
00	Premature infant	1.2%
0	Neonate	2.8%
1	Small Child	3.7%
1.5	Child	1.2%
2	Child	19.2%
3	Adult	37.5%
3.5	Adult, strong curve	1.8%
4	Large adult	22.6%
5	Extra-large adult	Minimal

Key Takeaways


- Blade sizing is related to patient size
- Not all configurations are available in all sizes

BLADE DESIGN & TECHNIQUE

Slight variations in laryngoscopic technique follow from the choice of blade design, and it is often a matter of personal preference.



Straight blade: Tip goes *under the epiglottis* and lifts it directly



Curved blade: Tip fits *into the vallecula* and indirectly lifts the epiglottis

VIDEO TIME OUT: 13 MIN

Anatomy & Direct Laryngoscopy (2:42 min)

<https://youtu.be/ThISkClbv7o>

Brief video that covers critical anatomy, blade types, and technique for direct laryngoscopy.



Direct and Indirect Laryngoscopy in 10 mins

<https://vimeo.com/129280563>


A clear, practical explanation of direct and indirect laryngoscopy, when each is preferred, pros/cons, trends in video use



HANDLE SHAPES




Volume by Size



Key Takeaways

- Almost all sales are either **Medium, Penlight, or Stubby**

LIGHTING SYSTEM: STANDARD VS. FIBEROPTIC LED



Key Takeaways

- Most disposables and reusables utilize an LED light source (bulb)
- The key difference is where the light source (bulb) is located
- Fiberoptic LED, often called Green or Greenline, always has a green collar
- Per ISO standards, Green Line blades and handles interchangeable regardless of manufacturer
- Although Fiberoptic LED is the most popular, Standard LED has superior light quality for fully-disposable products

VIDEO TIME OUT: 4 MIN STANDARD VS. FIBER OPTIC LED LIGHTING

Link here: <https://www.youtube.com/watch?v=isxV4pphogE>

- In under 4 mins, learn how to easily identify a Standard LED (aka Conventional) and Fiber Optic LED (Green Line)
- Have your IntuBrite Standard and Fiber Optic samples in hand to follow along



REUSABILITY


Reusable

Reusable items are designed to withstand reprocessing (high-level disinfection or sterilization). Most are guaranteed for duration (e.g. 5-7 years) and/or a defined number of autoclave cycles (e.g.) 4,000. Often, the electronic parts (bulb, battery, etc.) can be replaced when broken

Disposable

Disposable items are designed for single-patient use.

Intubations and Sales by Product Type



Key Takeaways

- All the money is in disposables.
- Winning disposable blades only to use with current reusable handles is a pot. wedge
- Winning reusables is unlikely. Don't devote time here.

REUSABLE VS DISPOSABLE: TRADE OFFS

Reusable Handle & Blade

Benefits

- High quality weight and feel
- What most physicians are accustomed to using – viewed as the gold standard

Drawbacks

- Performance failure/variability due to reprocessing wear & tear
 - Light quality degrades/ flicker
 - Battery failure, incorrect reinstallation
 - Component failure
- High initial capital investment
- Maintenance & replacement
- Reprocessing time & costs (chemicals, time, labor, supplies)

Disposable Handle & Blade

- Fresh battery and light for each use
- No maintenance or reprocessing resources
- Single patient use – no risk of cross-contamination

- Not perceived to have the quality of a reusable
- Cost is more 'tangible' to department
- Environmental impact perceived to be more than reusable

VIDEO REUSABLE & DISPOSABLE

Video laryngoscopes are also 'reusable' or 'disposable'

Reusable: Video monitor, handle and blade are intended for multi-patient use with handle cleaning and blade reprocessing between patients.

Disposable: These are really a combination of reusable and disposable components. Video monitor and camera are reusable. The 'blade' is either single use or has a disposable sheath.



**McGrath Mac
(Medtronic)**
Reusable monitor,
handle & metal stem
with camera.
Clear disposable
sheath fits over the
camera stem.



**King Vision
(Ambu)**
Reusable
monitor fits
into the
disposable
blade/handle
combo

OPTIONAL VIDEOS

Good views from practitioner's POV

https://www.youtube.com/watch?v=4V_poulbcnA

Direct Laryngoscopy <https://www.youtube.com/watch?v=AZeBumPaj4g>

Use of Mac and Miller blades for intubation

Anatomy and Intubation


<https://youtu.be/iOPpSGbuYmQ>

- Skip to the 6:30 min mark: @ 6:30 Anatomy; Great view of vallecula lift at 7:30 / Mallampati overview
- @ 8:30 Intubation Equipment; @ 12:00 Direct laryngoscopic intubation
The first 6.5 mins are pharmacology – double back to it if you can

Intubation using Direct and Indirect (Video) Laryngoscopes

<https://www.youtube.com/watch?v=gnkYGRMaw7o>

First 6 min shows use of different direct and video laryngoscopy devices



QUESTIONS?



 SunMed