

RESUS / HYPER / SAFET



# KEY PRODUCTS

- Resuscitation Devices
- Definitions
- VAP Overview
- AirFlow Manual Resuscitator

# DEFINITIONS

- **Manual Ventilation (also known as bag ventilation):**  
The intermittent manual compression of a gas-filled reservoir bag to force gases into a patient's lungs, this maintains oxygenation and carbon dioxide elimination during apnea or hypoventilation
- **Peak Inspiratory Pressure (PIP):**  
The highest level of pressure applied to the lungs during inhalation
- **Positive End-Expiratory Pressure (PEEP):**  
The pressure in the lungs that exists at the end of expiration
- **CPAP:**  
Continuous positive airway pressure
- **FiO<sub>2</sub> (Fraction of Inspired Oxygen):**  
The measurement of how much oxygen in the air is being inhaled from a ventilator
- **Tidal Volume:**  
The volume of air moved in or out of the lungs during quiet breathing




# PREVENTING VAP / VAE

## VENTILATION-ASSOCIATED PNEUMONIA (VAP) PREVENTION

ALSO KNOWN AS VENTILATION-ASSOCIATED EVENT (VAE)

### PATH TO PREVENTING VAP / VAE



### TYPICAL PATH TO VAP / VAE



# AIRFLOW MANUAL RESUSCITATORS

## Features

- Variety of styles
- Integrated Pressure Manometer
- Pliable body bag
- Highly customizable
- Bacterial/Viral Filter
- Single use, non-sterile, latex free

## Benefits

- Meet clinical preferences
- Apply appropriate PIP
- Protection for patients and clinicians





# RESUSCITATOR COMPONENTS

- Customize your bag!

## COMPONENTS AVAILABLE

Choose from the following components to create a custom resuscitation/ventilation bag

*(Some limitations apply)*



Integrated Pressure  
Manometer



Pop Off Valve



Bacterial / Viral  
Filter



CO<sub>2</sub> Indicator



PEEP Valve



Inflatable O<sub>2</sub> Bag



Corrugated  
O<sub>2</sub> Tubing



Collapsible O<sub>2</sub>  
Accumulator



Flex Neck




Hand Strap





Face Masks

# PRESSURE MANOMETER



# RESUSCITATOR COMPONENTS

- B/V Filters
- 19mm and 30mm options
- 99.99% efficiency
- Most popular – BF102

BACTERIAL / VIRAL FILTER			
CAT#	BF102	BF103	
DESCRIPTION	Accepts 19mm or 30mm	Accepts 19mm or 30mm	
B/V EFFICIENCY	>99.99%	>99.99%	
WEIGHT (G)	11.8	12.4	
ID / OD (MM)	19ID x 22ID / 30 OD	30ID x 22ID / 30OD	
DEAD SPACE (ML)	20	29	
RESISTANCE (CMH2O) 30LPM	<2.2	<2	
SAMPLING PORT	No	No	
PKG	50	50	



# RESUSCITATOR COMPONENTS



- Adjustable Peep Valves
- Pressure range: 5-20cmH2O
- Large print on scale for visibility
- Multiple options including MRI compatible

ITEM	DESCRIPTION	SIZE	PK
VP700	Adjustable PEEP Valve, 5-20 mg H2O	19 mm	20
VP702	Adjustable PEEP Valve + Adapter, Converts to 30 mm	19 mm + 30 mm Adapter	20
VP705	MRI Safe Adjustable PEEP Valve	19 mm	20

# RESUSCITATOR COMPONENTS

- CO2 Indicators

- **MaxCaps**

- For use up to 6 hours of continuous use with 360° view
- Efficient sponge material
- Neonatal, Pediatric, Regular sizing
- Fully compatible with SunMed devices
- Can be used with ET Tubes and AirQ



- **FloCap**

- For use up to 24 hours
- For patients greater than 15 kg
- Full compatible with all SunMed resuscitation devices
  - FloCap positioning is important to understand



# RESUSCITATOR COMPONENTS

- Masks

Vent Mask



Vent Mask II



Soft Seal



# RESOURCES

## RESUSCITATION / VENTILATION

### Ventlab AIRFLOW™ BAG - AF5140MB / AF5140MB-M5

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | CE MARK | DEHP



### INSTRUCTIONS FOR USE

#### PREPARATION FOR USE

1. Test the resuscitator for proper functioning: With patient port completely occluded, squeeze bag body to assure resistance is present. Positive needle movement should occur if equipped with a manometer.
2. Before using the mask, inspect for adequate inflation.
3. The manual resuscitator /ventilator may be used with a 19 mm peep valve with or without the expiratory filter attached. A 30 mm peep valve can be used with the expiratory filter attached. Attach peep valve to the expiratory filter that is attached to the exhalation port. Be sure that the accessory fits properly and does not interfere with compression of the resuscitator.
4. Actual PEEP may vary with patient lung compliance and resistance. Verify PEEP with a certified manometer.
5. For correct performance on the manual resuscitator/ventilator with oxygen reservoir, unfold the reservoir bag and assure that airflow is not restricted.

#### CAUTIONS

Expiratory filter, the filter must be dry and free of secretions. Wet filters have a high resistance that can impede ventilation and cause serious patient injury. Also, wet filters will not provide effective filtration.

#### DIRECTIONS FOR USE

1. Place the patient in a supine position. Establish and maintain an open airway.
2. Grasp the bag body with one hand.
3. Hold the mask between the index finger and thumb of the other hand. Place mask over face firmly to form a tight seal around the patient's nose and mouth.
4. Ventilate the patient by compressing the bag body for inhalation and releasing the bag body for patient's passive exhalation and bag body re-expansion. Continue this cycle as directed by medical authority.
5. Monitor peak airway pressures by observing the integrated manometer.  
**NOTE: If the patient is not intubated and patient condition allows, keep the peak inspiratory pressures <20 cm H<sub>2</sub>O which will help prevent gastric insufflation during mask ventilation.**
6. To remove vomitus: Disconnect resuscitator from patient. Tap the patient valve several times while squeezing the bag body. Re-test the resuscitator for proper functioning.
7. If patient is intubated, remove mask from patient port. Connect patient port directly to the endotracheal tube adapter. Continue ventilation.
8. To use supplemental oxygen: Connect oxygen supply tubing to O<sub>2</sub> source at appropriate flow rate. FDO<sub>2</sub> values may be affected if flow is not sufficient. Oxygen flow ≥ 15 lpm may be necessary. Do not let flow rate exceed 30 lpm due to possible increase in exhalation resistance.
9. When using the expiratory filter, monitor patient continuously while filter is in use. Please observe patient for proper chest movement during respiratory cycle. If ventilation is in question, remove filter from exhalation port and check filter for occlusion. If filter is occluded, discard and change filter.
10. Replace expiratory filter if used continuously for 24 hours, or more frequently, if resistance to flow reaches unacceptable levels.
11. Replace resuscitation bag when it is visibly soiled or per hospital policy — whichever comes first.

WARNING: Always read instructions for use on product packaging prior to use on of these devices.

- Resources can be found on the SunMed website [Sun-Med.com/resources](https://www.sunmed.com/resources)
- Product Literature
- Reference Material
- Mini-Brochures
- Videos
- And More!



Airway Management

Resuscitation/Ventilation

Anesthesia

Respiratory

Emergency

Oxygen Delivery

Diagnostics

Surgical

# KEY PRODUCTS

## Resuscitation Devices

- Hyperinflation System



# HYPERINFLATION COMPONENTS

## Features

- Integrated Pressure Manometer
- Pleated breathing bag
- Stay-put dial with red indication
- Highly customizable
- Swiveling patient port
- Single use, non-sterile, latex free

## Benefits

- Meet clinical preferences
- Apply appropriate static pressure
- Protection for patients



# HYPERINFLATION COMPONENTS

## SAFETY COMPONENTS



### INTEGRATED PRESSURE MANOMETER

Ventlab's innovative color-coded dial clearly denotes pressure delivery:

- Green** - Target pressure level <20 cm H<sub>2</sub>O
- Yellow** - Potentially dangerous pressure level
- Red** - Dangerous pressure level



### POP-OFF VALVES

40 cm H<sub>2</sub>O & 25 cm H<sub>2</sub>O  
Allows clinician to use an appropriate pop-off level

## CONVENIENCE COMPONENTS



### FACE MASKS

Inflatable and non-inflatable  
face masks available

- Customize your bag!

# RESOURCES

## HYPERINFLATION SYSTEM

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

The Hyperinflation System is available in five different pleated inflation bag sizes. Available with Ventlab's signature integrated pressure manometer and components to create a custom system.



- Resources can be found on the SunMed website [SunMed.com/resources](https://www.sunmed.com/resources)
- Product Literature
- Reference Material
- Mini-Brochures
- Videos
- And More!

### HYPERINFLATION SYSTEM

- System equipped with a stay-put dial to enable clinician to achieve targeted static pressure
- Five pleated breathing bag sizes available: 1/4, 1/2, 1, 2, 3 liter
- Seven inflatable face mask sizes available: small infant, infant, small child, child, adult (regular, large, wide)

ITEM#	DESCRIPTION	PK
HS1000	Hyperinflation System	20
HS4000	Hyperinflation System + Integrated Manometer	20

# KEY PRODUCTS

## NICU Products

- SafeT

# NICU OVERVIEW

- Neonatal Intensive Care Units have the most fragile patients in the entire hospital.
- These patients are rightfully protected by all, and SunMed products keeps their delicate anatomy top of mind.
- The SafeT is one of the safest ways to ventilate a neonate.
- The securement devices by Salter Labs keeps comfort and safety at the forefront for oxygen delivery.



# SAFET T-PIECE RESUSCITATOR



# RESUSCITATION REQUIREMENTS



# RESUSCITATION REQUIREMENTS

- During the first minute after delivery, the baby is being assessed
- If the baby is having difficulty breathing, persistent cyanosis, or HR < 100 bpm, then one of the following will be provided:
  - Positive Pressure Ventilation (PPV)
  - Continuous Positive Airway Pressure (CPAP)
- T-piece, BVM/Self-inflating bag, or flow-inflating bag will be used to provide the necessary support

# SAFET T-PIECE RESUSCITATOR

## Features


- Integrated Pressure Manometer
- PIP Controller
- Pop-off Valve
- Fits-All Oxygen Connector
- Expandable Circuit

## Benefits

- Meet clinical preferences
- Apply appropriate PIP and PEEP
- Adaptable in multiple scenarios



# SAFET T-PIECE RESUSCITATOR




Neonatal Resuscitator Features + Functionality	SafeT	BVM	Hyperinflation
Continuous O <sub>2</sub> Flow	●		●
Pre-Set + Deliver Consistent PIP	●		
Pre-Set + Deliver Consistent PEEP	●	●	●
Can Provide CPAP	●		●
Pre-set Maximum Pressure + Adjustability During Resus.	●		
Integrated Pressure Manometer	●	●	●
Pressure Relief Valve	●	●	●
Hand Compression Required		●	●
Portable	●	●	●



## SafeT™ T-PIECE RESUSCITATOR


LATEX FREE | SINGLE USE | DISPOSABLE

The SafeT T-Piece Resuscitator is the safe solution for neonatal resuscitation and CPAP. Delivering continuous oxygen flow while providing consistent PIP and PEEP pressures allows the clinicians to provide the necessary ventilatory support during resuscitation. SafeT helps prevent over-pressurization that can lead to barotrauma and pneumothorax in neonates.



### SafeT™ T-PIECE RESUSCITATOR

- Allows clinician to deliver consistent PIP and PEEP pressures
- Consistent PIP and PEEP ensures greater safety when ventilating babies' fragile lungs, particularly in an emergency situation
- Integrated pressure manometer provides clear visualization of ventilatory pressure ranges
- Pressure relief valve helps prevent over-pressurization of lungs
- Able to provide CPAP
- Includes a cushioned face mask for proper sealing during mask ventilation
- Easy to read controller and manometer makes set up simple and quick
- Universal red Fits-All O<sub>2</sub> connector allows immediate visualization of proper oxygen connection



CAT#	DESCRIPTION	MASK	PK
ST700	SafeT T-Piece Resuscitator	X-Small Neonate	10
ST701	SafeT T-Piece Resuscitator	Small Neonate	10
ST702	SafeT T-Piece Resuscitator	Small Infant	10
ST703	SafeT T-Piece Resuscitator	Infant	10
ST704	SafeT T-Piece Resuscitator	Neonate Round	10
ST705	SafeT T-Piece Resuscitator	Infant Round	10
ST707	SafeT T-Piece Resuscitator	-	10

- Resources can be found on the SunMed website [Sun-Med.com/resources](https://www.sun-med.com/resources)
- Product Literature
- Reference Material
- Mini-Brochures
- Videos
- And More!

# LIVE DEMONSTRATIONS

# AIRFLOW RESUS BAG



# HYPER RESUS BAG



# SAFET PIECE





# QUESTION AND ANSWER