



# Validating and comparing efficacy of supplemental oxygen delivery via four commercially available face masks

**Authors:** Herman Groepenhoff, Ph.D., Danielle Bayoro, M.D., Shiv Parikh, Dennis White, Edward A. Rose, M.D., Andreas Waldmann, MSc., and Michael J. Pedro, M.D., Vyair Medical, Mettawa, IL, USA

## Objectives


While closed simple oxygen masks deliver oxygen at high flow rates, the actual  $\text{FiO}_2$  received varies. Open oxygen masks allow increased comfort for the patient, improved communication and feeding, and higher clearance of carbon dioxide. This was an in vitro study of two open masks (Vyair AirLife™ Open, Southmedic OxyMask) and two closed masks (both Vyair Medical) at a variety of oxygen flow rates and tidal volumes in a simulated lung.

## Study methods

We tested three sizes each of four different types of masks. Masks were applied to dedicated head models that had internal tubing to simulate the nasal and oral cavities. The masks were attached to an adjustable oxygen supply. The nasal and oral tubes joined together and led to an ASL 5000 lung simulator and simulated tracheal  $\text{FiO}_2$  was continuously measured using an  $\text{O}_2$  sensor. Tests were repeated for oxygen flow rates between 1 and 15 LPM using simulated lung conditions. Each mask was tested nine times for 60 seconds per flow rate after 60 seconds washout.

## Results

At oxygen flow rates of 1 and 3 LPM, the open masks were able to deliver incremental  $\text{FiO}_2$  in excess of room air. At higher flow rates, the median  $\text{FiO}_2$  measured was similar across all masks, although the variability in measured  $\text{FiO}_2$  was significantly lower in the Open Mask than that seen in the closed masks ( $p < 0.001$ ). Variability in  $\text{FiO}_2$  was also significantly related to flow rate ( $p < 0.001$ ) in all masks but less prominently in the AirLife Open Mask. These patterns were seen in normal and obstructive lung settings.



Measured  $\text{FiO}_2$  by mask type at 15 LPM. Black horizontal lines represent median of the data. Each box represents the interquartile range (IQR). Upper and lower whiskers represent 1.5 times IQR. \*\*\*= $p < 0.001$

## What is an open mask?

Open oxygen masks have gained popularity due to their ability to deliver a consistent and reliable  $\text{FiO}_2$  at various flow rates with improved patient comfort, reduced rebreathing of carbon dioxide, and reduced hospital costs. Specifically, these devices can be used at a flow of 1–15 L/min to deliver 24–90% oxygen and were designed with one or more large openings to help clear carbon dioxide efficiently and prevent rebreathing. Additionally, the large holes were designed to improve activities of daily living by allowing the patient to drink and communicate.



Vyair AirLife™ Open Mask

## Take home message

The AirLife™ Open offers improved patient comfort and delivers oxygen at precise levels that are comparable to simple oxygen masks with decreased variability. The AirLife Open mask can be used at very low flow rates (<5 LPM) while the closed masks cannot and there is no risk of  $\text{CO}_2$  rebreathing with the open mask.

AirLife Open mask delivered significantly higher  $\text{FiO}_2$  than the Southmedic Oxymask. The open masks delivered  $\text{FiO}_2$  that was similar to that seen in the closed masks as flow rates increased beyond 5 LPM. There was a significantly higher variability of  $\text{FiO}_2$  in the closed masks compared to the open masks at higher oxygen flow rates.

### GLOBAL HEADQUARTERS

Vyair Medical, Inc. | 26125 North Riverwoods Blvd | Mettawa, IL 60045, USA

### For global distribution.

© 2021 Vyair. Vyair, the Vyair logo and all other trademarks or registered trademarks are property of Vyair Medical, Inc., or one of its affiliates. VYR-GBL-2100256