




enFlow™

Configuration guide

enFlow system has three main components:

1. Controller 980121EU
2. Warmer 980105VS
3. Disposable Parylene coated cartridge (with 3" extension) 980202EU or (without extension) 980200EU



Controller, PN 980121EU

The controller unit serves as the power supply for the warmer unit. It is designed to mount on an IV pole or sit on a table top. The front panel includes a temperature display and keypad.



Warmer holder, PN 980305VS, box of 20

The warmer holder affixes to the side of the controller to allow clinicians a place to hang the warmer when it is not in use.



Warmer, PN 980105VS

The warmer is designed to work in conjunction with the disposable cartridge to warm IV fluids. The innovative design of the enFlow warmer allows it to be placed within inches of the IV site, reducing the potential for fluid cooling within the IV line.



Cord clip, PN 980309VS-20, box of 20

The cord clip allows caregivers to affix the cord of the warmer to the patient's bed sheet or clothing.



Disposable cartridge, PN 980200EU, box of 30

The disposable cartridge can be connected to any standard luer IV set. The cartridge may stay in-line and travel with the patient for up to 24 hours. It requires less than 4 mL of priming volume.



Insulated warmer strap, PN 980304VS30, box of 30

The warmer strap with an integrated insulated pad allows the user to attach the warmer to the patient's limb. The insulated warmer strap should be used only when transporting the patient.



Disposable cartridge with 3" extension set, PN 980202EU, box of 30

Patient-dedicated cartridges with a 3"/7.5 cm extension set are also available for customers who require extra length at the end of the cartridge to allow for the placement of IV accessories.



Replacement pole clamp, PN 980330VS-1

A replacement pole clamp is available to secure the enFlow controller to an IV pole.



Power cords, P/N 91000170 (continental EU), 91000172 (UK)

The 10' power cord plugs into the back of the controller to supply power to the system.



System manual, PN 44000024

The system manual provides important information related to operation, service and preventative maintenance of the enFlow device. This is the English version. See below for additional languages.



IV clips, PN 980331VS-200, box of 200

IV line cord clips allow caregivers to affix the patient line to the warmer cord to prevent kinking in the line.








For ordering a ready to use enFlow system the following products need to be included








1. Controller 980121EU

PN – includes

- Controller unit – (sellable REF 980121EU)
- IV pole clamp – (spare part REF 980330VS-1)
- USB user manual

Power cord – Pick one:

Socket	Country	SKU
	Australia	AUS 91000180
	China	CN 91000179
	Brazil	BR 91000183
	Canada	USCAN - locking version 91000178USAL
	Japan	JP 91000182
	Continental Europe	Europe 91000170
	Denmark	DK 91000176

Socket	Country	SKU
	India	IN 91000175
	Israel	IL 91000174
	Italy	IT 91000173
	New Zealand	NZ 91000181
	UK	UK 91000172
	South Africa	RSA 91000177
	Switzerland	ZH 91000171

Paper Manual* – Pick one:

Language	SKU
English manual	44000024
German manual	44000073
Danish manual	44000074
Spanish manual	44000075
Finnish manual	44000076
French manual	44000077
Italian manual	44000078
Dutch manual	44000079
Norwegian manual	44000080
Swedish manual	44000081
Chinese manual	44000083
Turkish manual	44000084
Portuguese manual	44000085
Bulgarian manual	44000113
Croatian manual	44000114

Language	SKU
Czech manual	44000115
Greek manual	44000116
Polish manual	44000117
Russian manual	44000118
Serbian manual	44000119
Latvian manual	44000123
Lithuanian manual	44000124
Romanian manual	44000125
Slovakian manual	44000126
Slovenian manual	44000127
Hungarian manual	44000128
Estonian manual	44000140
Korean manual	44000141
Japanese manual	44000142
Traditional Chinese manual	44000143

2. Warmer

- 980105VS

3. Disposable parylene coated cartridge

- 980202EU – with 3" extension – order code includes 30 cartridges (EA)
- 980200EU – without extension – order code includes 30 cartridges (EA)

4. enCheck™

- enCheck, 980400
Recommended having one enCheck per customer or one per 10 enFlow systems. Used for annual preventive maintenance.

Optional products

- Warmer holder 980305VS, box of 20
Affixes to the side of the controller to store the warmer when not in use
- Warmer IV clip, 980331VS-200, box of 200
- Warmer cord clip, 980309VS-20, box of 20
- Insulated strap, 980304VS30, box of 30
- Controller pole clamp, 980330VS-1

GLOBAL HEADQUARTERS

Vyaire Medical, Inc.
26125 N. Riverwoods Blvd.
Mettawa, IL 60045
USA



Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands



Vyaire Medical, Inc.
26125 N. Riverwoods Blvd.
Mettawa, IL 60045
USA

AUSTRALIAN SPONSOR

Vyaire Medical Pty Ltd
Suite 5.03, Building C
11 Talavera Road
Macquarie Park, NSW 2113
Australia

*Mandatory in Europe
Performance is typical and depends upon end-to-end system configuration.
Actual performance may vary.

For EU, Canada and Australia distribution.

© 2021 Vyaire. Vyaire, the Vyaire logo and all other trademarks or registered trademarks are property of Vyaire Medical, Inc., or one of its affiliates. VYR-GBL-2100056

