

AIR-Q® AND AIR-Q®3 COMPARISON

The Air-Q®3 maintains the valuable features and offerings of the original Air-Q®, while introducing additional improvements.

AIR-Q®



Air-Q® Standard



Air-Q® Blocker



Air-Q®sp

AIR-Q®3 BENEFITS

- All silicone for ease of insertion, and softer mask is less irritating to the airway.
- Built up mask heel for improved seal.
- Dual gastric access for easy gastric suctioning.
- Larger bite block for improved patient comfort and airway protection.

AIR-Q®3



Air-Q®3



Air-Q®3G



Air-Q®sp3



Air-Q®sp3G

AIR-Q® TO AIR-Q®3 SIZE CROSS REFERENCES

AIR-Q SIZE	PATIENT WEIGHT (KGS)	MAX ETT (MM)	AIR-Q3 SIZE	PATIENT WEIGHT (KGS)	MAX ETT (MM)
-	-	-	0	< 2	3
-	-	-	0.5	2 - 4	4
1	< 7	4.5	1	4 - 7	4.5
1.5	7 - 17	5	1.5	7 - 17	5
2	17 - 30	5.5	2	17 - 30	5.5
2.5	30 - 50	6.5	3	30 - 60	7
3.5	50 - 70	7.5	4	60 - 80	8
4.5	70 - 100	8.5	5	> 80	9

AIR-Q® TO AIR-Q®3 CROSS REFERENCES

AIR-Q	AIR-Q3	AIR-Q BLOCKER	AIR-Q3G
10-3010	30105	10-5025	50305
10-3011	30105	10-5026	50305
10-3015	30155	10-5035	50405
10-3016	30155	10-5036	50405
10-3020	30205	10-5045	50505
10-3021	30205	10-5046	50505
10-3025	30305		
10-3026	30305		
10-3035	30405		
10-3036	30405		
10-3045	30505		
10-3046	30505		
AIR-QSP	AIR-QSP3	AIR-QSP	AIR-QSP3G
10-4010	40105	10-7025	60305
10-4015	40155	10-7026	60305
10-4020	40205	10-7035	60405
10-4025	40305	10-7036	60405
10-4035	40405	10-7045	60505
10-4045	40505	10-7046	60506
AIR-Q BLOCKER TUBE	AIR-Q3		
10-5027	No Cross		
10-5037	No Cross		
10-5047	No Cross		

AIR-Q® AND AIR-Q®3 REMOVAL STYLETS

ITEM	STYLET SIZE	AIR-Q3 SIZE	MAX. O.E.T.T. (MM)	WEIGHT (KGS)	PK
10-1006	#00 (small)	0.0 - 1.5	2.5 - 4.5	4-17	10
10-1005	#0 (medium)	2.0 - 3.0	4.5 - 6.0	17-60	10
10-1004	#1 (large)	4.0 - 5.0	6.0 - 9.0	> 60	10

FOR ADDITIONAL AIR-Q3 INFORMATION AND RESOURCES, VISIT: RESOURCES.SUN-MED.COM