


STYLETTES

LATEX FREE | SINGLE USE | DISPOSABLE

Stylettes are available in multiple sizes for both adult and pediatric patients. Choose from a variety of unique designs for a myriad of intubation applications.



ALUMINUM STYLETTE

- Malleable with rounded end
- Satin finish to reduce friction
- Pediatric stylette fits in 2-5 mm tubes
- Adult stylette fits in 6-11 mm tubes
- Individually wrapped

ITEM	DESCRIPTION	SIZE	PK
9-0204-11	Pediatric Stylette	11" x 6 Fr	25
9-0204-25	Adult Stylette	16" x 9 Fr	25

SUNGLIDE™ COATED STYLETTE

- Sterile
- Color coded adjustable stopper
- Coated surface reduces tension to decrease the need for lubricants

ITEM	DESCRIPTION	SIZE	PK
9-0209-71	SunGlide Stylette	13" x 6 Fr	25
9-0209-72	SunGlide Stylette	15.5" x 8 Fr	25
9-0209-73	SunGlide Stylette	18" x 10 Fr	25



LIGHTED STYLETTE

- Aids in blind intubation
- Adjustable stop

ITEM	DESCRIPTION	SIZE	PK
9-0210-01	Adult Stylette	15.75"	1

SUNSLIP™ CLEAR PVC COATED STYLETTE

- Sterile
- Malleable aluminum
- Sealed in a PVC sheath

ITEM	DESCRIPTION	SIZE	PK
9-0203-06	SunSlip with Sheath	10" x 6 Fr	20
9-0203-10	SunSlip with Sheath	11.5" x 10 Fr	20
9-0203-14	SunSlip with Sheath	12.5" x 14 Fr	20

STYLETTES

LATEX FREE

Stylettes are available in multiple sizes for both adult and pediatric patients. Choose from a variety of unique designs for a myriad of intubation applications.



STAINLESS STEEL STYLETTE

- Reusable
- Surgical stainless steel
- Malleable
- Non-sterile

ITEM	DESCRIPTION	SIZE	PK
9-0207-03	Adult Stylette	16.5" x 6 Fr	1
9-0207-01	Pediatric Stylette	11.5" x 6 Fr	1



FIBER OPTIC STYLETTE

- Reusable
- Adjustable stop
- Single use stops available

ITEM	DESCRIPTION	SIZE	PK
9-0204-18	Adult Stylette	18"	1
9-0204-14	Pediatric Stylette	14"	1
9-0240-59	Replacement Lamp	—	1
9-0210-05	Replacement Stops	—	5



SOLID COPPER STYLETTE

- Disposable
- Single Use
- Malleable
- Adjustable stop

ITEM	DESCRIPTION	SIZE	PK
9-0203-00	Solid Copper Stylette	14.75" x 8 Fr	1



HOLLOW COPPER STYLETTE

- Disposable
- Single Use
- Hollow copper tube features suction port

ITEM	DESCRIPTION	SIZE	PK
9-0206-01	Hollow Copper Stylette	16.5"	1