



UNI-LIM™ AND PEDI-LIM™ ANESTHESIA CIRCUIT

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

The single-limb circuit is the optimal choice for your Co-Axial Circle System needs. Designed to function as a standard circle anesthesia breathing system with streamlined tube-within-a-tube and eliminating bulky connectors.



UNI-LIM™ AND PEDI-LIM™ ANESTHESIA CIRCUITS

- Features a streamlined tube-within-a-tube design which eliminates twisting between the inspiratory and expiratory limbs
- Eliminates bulky connectors at the patient interface
- Provides natural heat and moisture conservation properties to that of a standard circle circuit
- Engineered with a secure-t inspiratory limb to prevent patient end disconnects – no floating internal tubes
- ISO-standard connections for secure, leak-resistant attachments

ITEM	DESCRIPTION	SIZE	ELBOW	FILTER	HME/HMEF	BAG	GSL	MASK
9061	60" Uni-Lim	Adult	GSE	(2) B/V	–	2 L	10' x .050" ID M/M	Large Adult 900 Series
8041	72" Uni-Lim	Adult	GSE	B/V	–	2 L	–	–
9432	96" Uni-Lim	Adult	GSE	–	Maxi HMEF	–	–	–
6754	60" Pedi-Lim	Pediatric	GSE	(2) B/V	–	2 L	10' x .050" ID M/M	–
8122	60" Pedi-Lim	Pediatric	GSE	B/V	–	1 L	–	–
9169	60" Pedi-Lim	Pediatric	GSE	B/V	–	2 L	10' x .050" ID M/M	–

CUSTOMIZE YOUR CIRCUIT

Choose from the following components to create a custom circuit to best meet your needs:



BAGS

Three standard bag size options available (1, 2, and 3 L) and a 2 liter EZ gripper bag that is easy to grip with a gloved hand



FILTERS

HMEs, HMEFs, bacterial viral, and HEPA filters with luer gas sampling ports available. Sampling ports located at the top of the housing, just as if a common gas sampling elbow were being used. ISO-standard connections for secure, leak-resistant attachment to breathing systems.



FLEX EXTENDERS

5.2" Flex Tubing with and without capped gas sampling port



FACE MASKS

Unique ergonomic mask shape that fits most facial configurations for improved sealing capabilities




GAS SAMPLING LINES

Standard and swivel lines available



WYES AND CONNECTORS

Elbows, parallel and swivel wyes available with and without port



YANKAUER

Sterile Yankauer with bulb tip and non-sterile Yankauer with regular tip available