



Premium face masks

Flex masks

Premium face masks

Our high-quality flexible face masks offer a soft face plate, anatomical teardrop shape and superior grip.

Cat. no.	Description Qty.	(per case)	Size	Weight (g)
6800EU	Clear, flexible dome, air-filled cushion face mask with no top valve	20	0	10.7–14.0
6810EU	Clear, flexible dome, air-filled cushion face mask with top valve	20	1	11.4–18.4
6820 6821	Clear, flexible dome, air-filled cushion face mask with top valve	20 100	2	12.9–22.3
6820BG 6821BG	Bubble gum scent, clear, flexible dome, air-filled cushion face mask with top valve	20 100	2	12.9–22.3
6830 6831	Clear, flexible dome, air-filled cushion face mask with top valve	20 100	3	28.0–37.0
6830BG 6831BG	Bubble gum scent, clear, flexible dome, air-filled cushion face mask with top valve	20 100	3	28.0–37.0
6840 6841	Clear, flexible dome, air-filled cushion face mask with top valve	20 100	4	39.0–56.8
6840BG 6841BG	Bubble gum scent, clear, flexible dome, air-filled cushion face mask with top valve	20 100	4	39.0–56.8
6850 6851	Clear, flexible dome, air-filled cushion face mask with top valve	20 100	5	45.0–67.7
6860 6861	Clear, flexible dome, air-filled cushion face mask with top valve	20 100	6	62.0–78.5

Materials	
Cushion	PVC Plastisol
Face plate	Polycarbonate
Scent additive	Bubble gum scent: 6820BG, 6830BG and 6840BG

Physical characteristics	
Color	Clear
Top valve	Size: 1, 2, 3, 4, 5 and 6
No valve	Size: 0

Connectors	
Face plate connector	ISO 15 mm male (ISO 5356-2) One on each size: 0, 1 and 2 ISO 22 mm female (ISO 5356-2) One on each size: 3, 4, 5 and 6
Prong rings	Attached to 22 mm female patient connection One on each: size: 3, 4, 5 and 6
Top valve	Luer connector (ISO 594)

Environmental conditions	
Operational	Temperature: -18 °C to 50 °C
Storage	Temperature: -40 °C to 60 °C

For more information, contact your local Vyaire representative.

GLOBAL HEADQUARTERS

Vyaire Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA



Vyaire Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

AUSTRALIAN SPONSOR

Vyaire Medical Pty Ltd
Suite 503, Level 5
11 Talavera Road
Macquarie Park, NSW, 2113
Australia

vyaire.com

For global distribution.

© 2022 Vyaire. Vyaire and the Vyaire Medical logo are trademarks or registered trademarks of Vyaire Medical, Inc or one of its affiliates. Medical devices class IIa according to Medical Devices Directive 93/42/EC. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labeling. VYR-GBL-2200149 1.0

