



Disposable face masks

Standard face masks

Standard disposable face masks

Cat. no.	Description	Size	Connector	Weight (g)	Qty. (cs.)
4810VS	Anesthesia face mask, neonatal, air filled cushion, Red hook ring, with top valve	1	15 mm	8	20
4820VS	Anesthesia face mask, infant, air filled cushion, Green hook ring, with top valve	2	15 mm	11	20
4830VS	Anaesthesia face mask, toddler, air filled cushion, Purple hook ring, with top valve	3	22 mm	16	20
4840VS	Anesthesia face mask, child, air filled cushion, Yellow hook ring, with top valve	4	22 mm	20	20
4852VS	Anesthesia face mask, small adult, air filled cushion, White hook ring, with top valve	5	22 mm	27	20
4860VS	Anesthesia face mask, Medium adult, air filled cushion, Blue hook ring, with top valve	6	22 mm	29	20
4870VS	Anesthesia face mask, Large adult, air filled cushion, Orange hook ring, with top valve	7	22 mm	36	20

Materials

Cushion	Polyvinyl chloride (PVC)
Face plate	Polycarbonate
Hook ring	Polypropylene
Valve	Silicone
Top valve	Luer connector

Physical characteristics

Mask color	Blow molded clear cushion and face plate
Top valve	All sizes have a top valve with blue silicone stem

Key benefits:

- Semi-flexible mask
- Anatomical
- Top valve
- Colored hook ring for identification
- Free from natural rubber latex and DEHP

For more information, contact your local Vyaire representative.

GLOBAL HEADQUARTERS

Vyaire Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA



Vyaire Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

AUSTRALIAN SPONSOR

Vyaire Medical Pty Ltd
Suite 503, Level 5
11 Talavera Road
Macquarie Park, NSW, 2113
Australia

vyaire.com

For global distribution.

© 2022 Vyaire. Vyaire and the Vyaire Medical logo are trademarks or registered trademarks of Vyaire Medical, Inc or one of its affiliates. Medical devices class IIa according to Medical Devices Directive 93/42/EC. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labeling. VYR-GBL-2200150 1.0

