



Vyaire™ oxygen therapy consumables

Product portfolio

Vyaire™ oxygen therapy consumables product portfolio

Our extensive portfolio of oxygen therapy solutions empowers you to provide the highest quality patient care. From the most basic nasal cannula to a 3 in 1 non-rebreather mask, Vyaire products meet the needs of oxygen care patients both in the hospital and at home.

Why Vyaire oxygen therapy solutions?

- Clinical differentiation through an industry-leading and unique product portfolio
- Standardisation around a core respiratory portfolio with a leading set of products
- Quality care solutions that reduce variability in processes and products

Table of contents

Oxygen cannulas	3–4
Soft cushion cannulas.....	3
Foam cover cushion cannulas	3
Cushion cannulas	4
Standard nasal cannulas	4
Oxygen masks	5–7
Disposable O ₂ masks	5
Venturi masks	6
Oxygen & Aerosol masks	7
Tracheostomy masks	7
Bubble humidifiers and large volume nebulizer	8–9
Bubble humidifier	8
Misty Ox™	8
Oxygen tubing	9

Oxygen cannulas

All oxygen cannulas come with crush-resistant tubing to minimize occlusion and kinking. Select cannulas come with U/Connect-It® oxygen tubing for easy connection to any oxygen source.

Soft cushion cannulas

Our soft cushion cannulas include lariat tubing that is approximately 37% softer than our standard oxygen cannulas.* This soft tubing helps reduce chafing behind the ear to minimise skin irritation. They are made without DEHP.

Soft cushion cannula		
Catalog no	Description	Qty. (cs.)
SFT2611	Neonate, 2.1 m (7') O ₂ tubing	50
SFT2601	Infant, 2.1 m (7') O ₂ tubing	50
SFT2602	Paediatric, 2.1 m (7') O ₂ tubing	50
SFT2692	Paediatric, 2.1 m (7') O ₂ U/Connect-It® tubing	50
SFT2604	Adult, 1.2 m (4') tubing	50
SFT2699	Adult, 2.1 m (7') O ₂ U/Connect-It tubing	50
SFT2614U	Adult, 4.3 m (14') O ₂ U/Connect-It tubing	50
SFT2625	Adult, 7.8 m (25') O ₂ tubing	25



Soft cushion cannula

Foam cover cushion cannulas

Our cushion cannulas with soft foam ear covers help reduce skin irritation around the ears. The soft, gray covers slide smoothly on the headset tubing, making them easy to position for enhanced comfort and fit. They are made without DEHP.

Foam cover cushion cannula		
Catalog no	Description	Qty. (cs.)
FM2692	Paediatric, 2.1 m (7') O ₂ U/Connect-It tubing	50
FM2606	Adult, 2.1 m (7') no tubing	50
FM2604	Adult, 1.2 m (4') tubing	50
FM2699	Adult, 2.1 m (7') O ₂ U/Connect-It tubing	50
FM2614	Adult, 4.3 m (14') O ₂ tubing	50
FM2625	Adult, 7.8 m (25') O ₂ tubing	25
002016	Cannula ear covers (by pair)	50
SFT2614U	Adult, 4.3 m (14') O ₂ U/Connect-It tubing	50
SFT2625	Adult, 7.8 m (25') O ₂ tubing	25



Foam cover cushion cannula


* Based on a material durometer difference of 30 using the Shore A hardness scale

Oxygen cannulas (cont'd)

Cushion cannulas

Our cushion cannulas have a soft, clear cushion design that supports the nasal prongs to optimise patient comfort. Cushion cannulas are available in sizes from neonate to adult with a variety of oxygen tubing lengths. They are made without DEHP.

Cushion cannula		
Catalog no	Description	Qty. (cs.)
002611	Neonate, 2.1 m (7') O ₂ tubing	50
002601	Infant, 2.1 m (7') O ₂ tubing	50
002662	Paediatric, 2.1 m (7') no tubing	50
002602	Paediatric, 2.1 m (7') O ₂ tubing	50
002692	Paediatric, 2.1 m (7') O ₂ U/Connect-It tubing	50
002699	Adult, 2.1 m (7') O ₂ U/Connect-It tubing	50
002600-14U	Adult, 4.3 m (14') O ₂ U/Connect-It tubing	50
002600-25	Adult, 7.8 m (25') O ₂ tubing	25
002600-50	Adult, 15 m (50') O ₂ tubing	20
002600-4	Adult, 1.2 m (4') tubing	50
002630-4	Adult nasal cannula, 1.2 m (4') Tubing	50



Cushion cannula

Standard nasal cannulas

Our standard nasal cannulas offer an angulated, flexible lip plate, which provides comfort at the tip of the nasal region with beveled tips to provide a softer edge. Disposable cannulas are curved and available with flared or non-flared tips. They are packaged separately with varying lengths of oxygen tubing. They are made without DEHP.

Standard cannula		
Catalog no	Description	Qty. (cs.)
001365	Adult, flared tip, 2.1 m (7') O ₂ U/Connect-It tubing	50
001321	Adult, flared tip, 4.3 m (14') O ₂ tubing	50
001322	Adult, flared tip, 6.4 m (21') O ₂ tubing	25
001366	Adult, non-flared tip, 2.1 m (7') O ₂ U/Connect-It tubing	50
001327	Adult, non-flared tip, 6.4 m (21') O ₂ tubing	25
001329	Adult, non-flared tip, 15 m (50') O ₂ tubing	15



Standard cannula

Oxygen masks

Vyaire oxygen masks are latex-free and made from soft, high-grade vinyl resins. The smooth, feathered edges are designed for patient comfort while reducing irritation points. Our oxygen masks are designed for anatomical compatibility. A full range is available to meet all requirements in infant, paediatric and adult sizes. Select cannulas come with U/Connect-It oxygen tubing for easy connection to any oxygen source. The three-in-one masks are packaged as high-concentration, non-rebreathing masks complete with reservoir bags and one-way valves.



Oxygen mask		
Catalog no	Description	Qty. (cs.)
001262	Paediatric, medium-concentration, 2.1 m (7') O ₂ tubing	50
001267	Paediatric, high-concentration, rebreather, 2.1 m (7') crush-resistant O ₂ tubing	50
001268	Paediatric, high-concentration, non-rebreather 3-in-1 with safety vent, 2.1 m (7') crush-resistant O ₂ tubing	50
001269	Paediatric, high-concentration, non-rebreather, 2.1 m (7') O ₂ tubing	50
001260	Paediatric, medium-concentration, short style, 2.1 m (7') O ₂ tubing	50
001264	Paediatric short style oxygen mask, high concentration rebreather, mask only	50
001265	Paediatric, high-concentration, short style, 2.1 m (7') crush-resistant O ₂ tubing	50
001200	Adult oxygen mask, under-the-chin style, medium concentration mask only	50
001201	Adult, medium-concentration, 2.1 m (7') crush-resistant O ₂ tubing	50
001207	Adult high-concentration mask, under-the-chin style, mask only	50
001202	Adult oxygen mask, 3 in 1 mask with safety vent, mask only	50
001204	Adult oxygen mask, high concentration rebreather mask (no valves), mask only	50
001361	Adult, medium-concentration, 2.1 m (7') crush-resistant O ₂ U/Connect-It tubing	50
001203	Adult, high-concentration, non-rebreather 3-in-1 with safety vent, 2.1 m (7') crush-resistant O ₂ tubing	50
001362	Adult, high-concentration, non-rebreather 3-in-1 with safety vent, 2.1 m (7') crush-resistant O ₂ U/Connect-it tubing	50
001210	Adult, high-concentration, non-rebreather 3-in-1, 7' crush-resistant O ₂ tubing	50
001211	Adult, high-concentration, non-rebreather 3-in-1 with safety vent, 7' O ₂ tubing	50
001212	Adult, high-concentration, non-rebreather 3-in-1, 2.1 m (7') O ₂ tubing	50
001205	Adult oxygen, high concentration rebreather mask, 2.1 m (7') crush resistant oxygen tubing	50
001220	Face tent mask vinyl	50

Oxygen masks (cont'd)

Venturi masks

We offer a wide selection of both standard and Dual Dial™ Venturi-style masks that feature either six diluter jets or two oxygen concentration dials. All masks are under-the-chin style and supplied with a six-inch flextube, humidification hood and seven-foot crush-resistant tubing. We also offer a Venturi-style with a tracheostomy mask for use with a 50 psi oxygen source and flowmeter.



001240 Adult Venturi-style mask

Venturi mask		
Catalog no	Description	Qty. (cs.)
001250	Paediatric, 15 cm (6") flextube, 6 diluter jets (24%, 28%, 31%, 35%, 40%, 50%), humidification hood, 2.1 m (7') crush-resistant O ₂ tubing	50
001240	Adult, 15 cm (6") flextube, 6 diluter jets (24%, 28%, 31%, 35%, 40%, 50%), humidification hood, 2.1 m (7') crush-resistant O ₂ tubing	50
001248	Adult, tracheostomy, 15 cm (6") flextube, 6 diluter jets (24%, 28%, 31%, 35%, 40%, 50%), humidification hood, 2.1 m (7') crush-resistant O ₂ tubing	50
001363	Adult, 15 cm (6") flextube, 6 diluter jets (24%, 28%, 31%, 35%, 40%, 50%), humidification hood, 2.1 m (7') crush-resistant O ₂ U/Connect-It tubing	50
001241	Adult Venturi mask, 15 cm (6") flextube, 24% oxygen diluter jet, humidification hood and 2.1 m (7') crush-resistant tubing	50
001255	Adult, 15 cm (6") flextube, Dual Dial, humidification hood, 2.1 m (7') crush-resistant O ₂ tubing	50
001247	Venturi tracheostomy tee with 15 cm (6") flextube, 6 diluter jets (24%, 28%, 31%, 35%, 40% and 50% oxygen concentrations), humidification hood, 2.1 m (7') crush-resistant tubing	50

Oxygen masks (cont'd)

Oxygen & Aerosol masks

These soft, high-quality disposable aerosol masks are made from high-grade resins. The smooth, feathered edges are designed for patient comfort by reducing irritation points. All our masks are designed for anatomical compatibility. A comprehensive line is available to meet requirements in paediatric and adult sizes.

Aerosol mask		
Catalog no	Description	Qty. (cs.)
001263	Paediatric	50
001264	Paediatric short style oxygen mask, high concentration rebreather, mask only	50
001206	Adult	50
001209	Adult aerosol mask, 1.5 m (5') corrugated tubing	50
001222	Adult face tent mask, 1.5 m (5') corrugated tubing	50
001221	Adult face tent mask, 1.8 m (6') corrugated tubing	50



001206 Adult aerosol mask

Tracheostomy masks

The disposable tracheostomy masks are available in paediatric and adult sizes, and incorporate several patient features. The neckband is made from a comfortable, nonbiting material, and the swivel tubing connector permits access from either side of the patient. The special strap clips allow the mask to be removed with minimal disturbance to the patient.

Tracheostomy masks		
Catalog no	Description	Qty. (cs.)
001225	Paediatric	50
001226	Adult	50
001227	Adult, large	50




001226 Adult tracheostomy mask

Oxygen tubing

Our extensive line of disposable vinyl-tipped oxygen tubing comes with crush-resistant lumen to resist occlusion. Choose from a variety of lengths ranging from 2.1 m (7") to 15 m (50'). Tubing 001350 also comes with a U/Connect-It connector to attach directly to flowmeters without an adapter.

Vyaire disposable corrugated tubing for aerosol and IPPB applications offers flexibility and convenience to meet therapeutic needs. Cuffs on all tube ends lock to 22 mm (7/8") connections.

Oxygen tubing		
Catalog no	Description	Qty. (cs.)
001302*	2.1 m (7') crush resistant	50
001350*	2.1 m (7') crush resistant with U/Connect-It	50
001303*	4.3 m (14') crush resistant	50
001304	Airlife® oxygen tubing 6.4 m (21') crush resistant	25
001305*	7.6 m (25') crush resistant	25
001307	10 m (30') crush resistant	20
001308*	12 m (40') crush resistant	20
001306*	15 m (50') crush resistant	15
001303GRN	4.3 m (14') crush resistant, green	50
001304GRN	Tubing 6.4 m (21') crush resistant green	25
001308GRN	12 m (40') crush resistant green	20
001305GRN	27.6 m (25') crush resistant, green	25
001306GRN	15 m (50') crush resistant, green	15
001330	Tubing oxygen 2.1 m (7') smooth bore	50



Oxygen tubing

Corrugated tubing

Catalog no	Description	Qty. (cs.)
001400*	1.8 m (6')	50
001410*	1.8 m (6') in 6" segments	50
001404*	30.5 m (100') in 6" segments	1
001405*	30.5 m (100') in 6" segments with dispenser box	1
001421	5' blue in 6" segments	50
001420	1.8 m (6') blue in 6" segments	50
001426	30.5 m (100') blue in 6" segments	1
001427	30.5 m (100') blue in 6" segments with dispenser box	1
001450	Corrugated flexible polyethylene and EVA tubing, segmented, 1.8 m (6') length, 22 mm I.D.	

* Does not contain DEHP.



Corrugated tubing

**Products may not be available in all regions.
For more information, contact your local Vyaire
representative.**

GLOBAL HEADQUARTERS

Vyaire Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA



Vyaire Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA



Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands

AUSTRALIAN SPONSOR

Vyaire Medical Pty Ltd
Suite 503, Level 5
11 Talavera Road
Macquarie Park, NSW, 2113
Australia

vyaire.com

For European Union distribution only.

Trademarks are the property of their respective owners.

© 2022 Vyaire. Vyaire, the Vyaire logo and AirLife are trademarks or registered trademarks of Vyaire Medical, Inc., or one of its affiliates. Medical devices class IIa according to Medical Devices Directive 93/42/EC.

Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-GBL-2200159

