



Diagnostic Cardiology, IEC

Validated accessories for GE systems

Validated accessories for GE diagnostic cardiology

Our Diagnostic Cardiology cables, leadwires, connectors and papers are official, validated supplies for your GE Healthcare equipment. They have been designed for optimal performance and use with GE Healthcare ECG systems, connecting the device to the patient to capture an accurate reading.

These validated accessories include:

- Cables and Leadwires
- Connectors and adapters
- Paper
- Other accessories such as gels and electrodes

For the following systems:

- Cardiosoft CASE
- MAC600
- MAC800
- MAC1200
- MAC2000
- MAC3500
- MAC5500 HD
- SEER Light
- SEER 1000
- KISS Suction Pump
- Cardiomem® CM 3000 SM, SM/A, Z, ZSM
- Cardiomem® CM 3000-12(BT)
- SEER® 12

Key Features

- Validated for use with GE Diagnostic Cardiology Systems
- Shielded to help reduce electromagnetic interference
- Designed to securely connect to electrodes to prevent accidental disconnection



For more information or to order product, contact:

International Customer Support

International Customer Support: +44-1256-588062
Email: INTL-CustomerService@vyaire.com
(US) 1-833-327-3284, email: customersupport@vyaire.com

UK Ordering and Customer Support

UK Freephone Customer Support: 0800 051 7201 (Option 2)
Email: GMB-UK-ARD-Customer-Service@vyaire.com



Device overview



2016560-001, cable, CAM-14



22341808V, cable, MAC



2009828-024, paper, DCAR



900178-101, connector, banana




900178-202, connector, alligator




900178-103, connector, grabber



2066867-014, adapter clip, wide mouth



2089797-003, leadwire set, universal




2016033-002, leadwire set, banana

Table of contents

CardioSoft & CASE	5	MAC 1600	10	SEER Light	15
ECG trunk cables.....	5	ECG trunk cables.....	10	ECG leadwire sets	15
ECG leadwire sets.....	5	ECG trunk cables with integrated leadwires.....	10	Accessories.....	15
Individual ECG leadwires.....	5	ECG leadwire sets	10	SEER 1000	16
Banana connectors.....	5	Banana connectors.....	10	ECG leadwire sets	16
Universal adapter clips.....	5	Universal adapter clips.....	10	KISS system	17
Mactrode alligator clips.....	5	Mactrode alligator clips.....	10	ECG cables for KISS systems.....	17
Grabber adapter clips.....	5	Grabber adapter clips.....	10	ECG leadwires for KISS system	17
ECG accessories	6	ECG accessories	10	KISS leadwires accessories.....	17
ECG paper	6	ECG paper	10	Accessories.....	17
MAC 600		MAC 2000	11	Cardiomem® CM 3000 SM, SM/A, Z, ZSM; Cardiomem® CM 3000- 12(BT); SEER(R) 12	18
ECG trunk cables.....	7	ECG trunk cables.....	11	5-lead combined cable & leadwire...18	
ECG trunk cables with integrated leadwires.....	7	ECG trunk cables with integrated leadwires.....	11	5-lead Separate cable & leadwire18	
ECG leadwire sets.....	7	ECG leadwire sets	11	7-lead Combined cable & leadwire..18	
Banana connectors.....	7	Banana connectors.....	11	7-lead Separate cable & leadwire.....18	
Universal adapter clips.....	7	Universal adapter clips.....	11	12-lead for CM 3000-12(BT) and SEER 12.....	18
Mactrode alligator clips.....	7	Mactrode alligator clips.....	11	Electrodes	19
Grabber adapter clips.....	7	Grabber adapter clips.....	11		
ECG accessories	7	ECG accessories	11		
ECG paper	7	ECG paper	11		
MAC 800	8	MAC 3500	12		
ECG trunk cables.....	8	ECG trunk cables.....	12		
ECG trunk cables with integrated leadwires.....	8	ECG leadwire sets	12		
ECG leadwire sets.....	8	Banana connectors.....	12		
Banana connectors.....	8	Universal adapter clips.....	12		
Universal adapter clips.....	8	Mactrode alligator clips.....	12		
Mactrode alligator clips.....	8	Grabber adapter clips.....	12		
Grabber adapter clips.....	8	ECG accessories	12		
ECG accessories	8	ECG paper	12		
ECG paper	8	MAC 5500 HD	13		
MAC 1200	9	ECG trunk cables.....	13		
ECG trunk cables.....	9	ECG leadwire sets	13		
ECG trunk cables with integrated leadwires.....	9	Individual ECG leadwires.....	13		
ECG leadwire sets.....	9	Banana connectors.....	13		
Banana connectors.....	9	Universal adapter clips.....	13		
Universal adapter clips.....	9	Mactrode alligator clips.....	13		
Mactrode alligator clips.....	9	Grabber adapter clips.....	13		
Grabber adapter clips.....	9	ECG accessories	14		
ECG accessories	9	ECG paper	14		
ECG paper	9				

CardioSoft® & CASE® Stress Test Systems

Part number	Product description	Order UOM (unit of measure)	Visual aid
ECG trunk cables			
2016560-002	Cable, DCAR, CAM-14 Cable for use with CardioSoft and CASE, 4.6 m	Each	
ECG leadwire sets			
2089797-004	Leadwire set, DCAR universal, mixed length for 12-lead. Set of 10 leads (6 leads 66 cm, 2 leads 90 cm, 2 leads 100 cm long)	Each (set of 10 mixed lengths)	
2089797-003	Leadwire set, DCAR universal, mixed length for 15-lead. Set of 14 leads (10 leads 66 cm, 2 leads 90 cm, 2 leads 100 cm long)	Each (set of 14 mixed lengths)	
38401816	Leadwire, DCAR, banana-ended, 12-lead. Set of 10 leads (4 leads 1295 mm long, 6 leads 737 mm long)	Each (set of 10 mixed lengths)	
2016033-002	Leadwire, DCAR, banana-ended, for 15-lead, mixed lengths. Set of 14 leads (2 leads 660 mm long, 2 leads 914 mm long, 10 leads 1016 mm long)	Each (set of 14 mixed lengths)	
Individual ECG leadwires			
2089799-004	Leadwire, DCAR universal, C3R, 66 cm	Each	
2089799-005	Leadwire, DCAR universal, C4R, 66 cm	Each	
2089799-006	Leadwire, DCAR universal, C7, 66 cm	Each	
Banana connectors for use with universal leadwires			
900178-101	Connector, DCAR, banana for 12-lead	Each (set of 10)	
900178-201	Connector, DCAR, banana for 15-lead	Each (set of 14)	
Universal adapter clips for use with banana-ends			
2056813-010	Adapter clip, DCAR, Sure-Lock clip for 12-lead	Package (10 each)	
2056813-014	Adapter clip, DCAR, Sure-Lock clip for 15-lead	Package (14 each)	
2066867-010	Adapter clip, DCAR, wide mouth clip for 12-lead	Package (10 each)	
2066867-014	Adapter clip, DCAR, wide mouth clip for 15-lead	Package (14 each)	
Mactrode alligator clips for use with banana-ends			
900178-102	Adapter clip, DCAR, Mactrode alligator clip for 12-lead	Each (set of 10)	
900178-202	Adapter clip, DCAR, Mactrode alligator clip for 15-lead	Each (set of 14)	
Grabber adapter clips for use with banana-ends			
900178-103	Connector, DCAR, grabber for 12-lead	Each (set of 10)	
900178-203	Connector, DCAR, grabber for 15-lead	Each (set of 14)	

Part number	Product description	Order UOM <i>(unit of measure)</i>	Visual aid
ECG accessories			
2001926-001	Accessory, DCAR, leadwire blank plug for CAM-14	Pack (set of 4)	
401089-003	Accessory, DCAR, leadwire separators, set of 6	Each (set of 6)	
ECG paper for use with CASE® system only			
2009828-021	ECG paper, A4, red grid, 10 year archival life, z-fold	Case (2400 sheets) ((8 packages of 300 sheets))	
2009828-020	ECG paper, letter, red grid, 10 year archival life, z-fold	Case (2400 sheets) ((8 packages of 300 sheets))	
2009828-034	ECG paper, letter, red grid, 25 year archival life, Z-fold	Case (2400 sheets) ((8 packages of 300 sheets))	
2009828-035	ECG paper, A4, red grid, 25 year archival life, z-fold	Case (2400 sheets) ((8 packages of 300 sheets))	

MAC 600

Part number	Product description	Case qty	Visual aid
ECG trunk cables			
22341808V	Cable, DCAR, 12-lead for MAC systems, 210 cm	Each	
ECG trunk cables with integrated leadwires			
2029890-001V	Cable trunk MAC500/1200 10-lead	Each	
ECG leadwire sets			
2089797-004	Leadwire set, DCAR universal, mixed length for 12-lead. Set of 10 leads (6 leads 66 cm, 2 leads 90 cm, 2 leads 100 cm long)	Each (set of 10)	
38401816	Leadwire, DCAR, banana-ended, 12-lead Set of 10 leads (4 leads with length of 1295 mm, 6 leads with length of 737 mm)	Each (set of 10)	
Banana connectors for use with universal leadwires			
900178-101	Connector, DCAR, banana for 12-lead	Each (set of 10)	
Universal adapter clips for use with banana-ends			
2056813-010	Adapter clip, DCAR, Sure-Lock clip for 12-lead	Package (10 each)	
2066867-010	Adapter clip, DCAR, wide mouth clip for 12-lead	Package (10 each)	
Mactrode alligator clips for use with universal leadwires			
900178-102	Adapter clip, DCAR, Mactrode alligator clip for 12-lead	Each (set of 10)	
Grabber adapter clips for use with universal leadwires			
900178-103	Connector, DCAR, grabber for 12-lead	Each (set of 10)	
ECG accessories			
401089-003	Accessory, DCAR, leadwire separators	Each (set of 6)	
ECG paper			
2030887-001	Paper, DCAR, full red grid, Z-fold, 80 mm	Case (2800 sheets) (10 packages of 280 sheets)	

MAC 800

Part number	Product description	Case qty	Visual aid
ECG trunk cables			
22341808V	Cable, DCAR, 12-lead for MAC systems, 210 cm	Each	
ECG trunk cables with integrated leadwires			
2029890-001V	Cable trunk MAC500/1200 10-lead	Each	
ECG leadwire sets			
2089797-004	Leadwire set, DCAR universal, mixed length for 12-lead. Set of 10 leads (6 leads 66 cm, 2 leads 90 cm, 2 leads 100 cm long)	Each (set of 10)	
38401816	Leadwire, DCAR, banana-ended, 12-lead Set of 10 leads (4 leads with length of 1295 mm, 6 leads with length of 737 mm)	Each (set of 10)	
Banana connectors for use with universal leadwires			
900178-101	Connector, DCAR, banana for 12-lead	Each (set of 10)	
Universal adapter clips for use with banana-ends			
2056813-010	Adapter clip, DCAR, Sure-Lock clip for 12-lead	Package (10 each)	
2066867-010	Adapter clip, DCAR, wide mouth clip for 12-lead	Package (10 each)	
Mactrode alligator clips for use with universal leadwires			
900178-102	Adapter clip, DCAR, Mactrode alligator clip for 12-lead	Each (set of 10)	
Grabber adapter clips for use with universal leadwires			
900178-103	Connector, DCAR, grabber for 12-lead	Each (set of 10)	
ECG accessories			
401089-003	Accessory, DCAR, leadwire separators	Each (set of 6)	
ECG paper			
2036970-001	ECG paper, red grid, 8 year archival life, z-fold	Case (2000 sheets) (10 packages of 200 sheets)	

MAC 1200

Part number	Product description	Case qty	Visual aid
ECG trunk cables			
22341808V	Cable, DCAR, 12-lead for MAC systems, 210 cm	Each	
ECG trunk cables with integrated leadwires			
2029890-001V	Cable trunk MAC500/1200 10-lead	Each	
ECG leadwire sets			
2089797-004	Leadwire set, DCAR universal, mixed length for 12-lead. Set of 10 leads (6 leads 66 cm, 2 leads 90 cm, 2 leads 100 cm long)	Each (set of 10)	
38401816	Leadwire, DCAR, banana-ended, 12-lead Set of 10 leads (4 leads with length of 1295 mm, 6 leads with length of 737 mm)	Each (set of 10)	
Banana connectors for use with universal leadwires			
900178-101	Connector, DCAR, banana for 12-lead	Each (set of 10)	
Universal adapter clips for use with banana-ends			
2056813-010	Adapter clip, DCAR, Sure-Lock clip for 12-lead	Package (10 each)	
2066867-010	Adapter clip, DCAR, wide mouth clip for 12-lead	Package (10 each)	
Mactrode alligator clips for use with universal leadwires			
900178-102	Adapter clip, DCAR, Mactrode alligator clip for 12-lead	Each (set of 10)	
Grabber adapter clips for use with universal leadwires			
900178-103	Connector, DCAR, grabber for 12-lead	Each (set of 10)	
ECG accessories			
401089-003	Accessory, DCAR, leadwire separators	Each (set of 6)	
ECG paper			
2009828-061	ECG paper, letter, red grid with blank data area, 10 year archival life, z-fold	Case (2400 sheets) (16 packages of 150 sheets))	
22616605	ECG paper, A4, red grid, 8 year archival life, z-fold	Case (1500 sheets) (10 packages of 150 sheets))	
22616611	ECG paper, A4, red grid, 25 year archival life, z-fold		

MAC 1600

Part number	Product description	Case qty	Visual aid
ECG trunk cables			
22341808V	Cable, DCAR, 12-lead for MAC systems, 210 cm	Each	
ECG trunk cables with integrated leadwires			
2029890-001V	Cable trunk MAC500/1200 10-lead	Each	
ECG leadwire sets			
2089797-004	Leadwire set, DCAR universal, mixed length for 12-lead. Set of 10 leads (6 leads 66 cm, 2 leads 90 cm, 2 leads 100 cm long)	Each (set of 10)	
38401816	Leadwire, DCAR, banana-ended, 12-lead Set of 10 leads (4 leads with length of 1295 mm, 6 leads with length of 737 mm)	Each (set of 10)	
Banana connectors for use with universal leadwires			
900178-101	Connector, DCAR, banana for 12-lead	Each (set of 10)	
Universal adapter clips for use with banana-ends			
2056813-010	Adapter clip, DCAR, Sure-Lock clip for 12-lead	Package (10 each)	
2066867-010	Adapter clip, DCAR, wide mouth clip for 12-lead	Package (10 each)	
Mactrode alligator clips for use with universal leadwires			
900178-102	Adapter clip, DCAR, Mactrode alligator clip for 12-lead	Each (set of 10)	
Grabber adapter clips for use with universal leadwires			
900178-103	Connector, DCAR, grabber for 12-lead	Each (set of 10)	
ECG accessories			
401089-003	Accessory, DCAR, leadwire separators	Each (set of 6)	
ECG paper			
2009828-060	ECG paper, letter, red grid, 10 year archival life, z-fold	Case (2400 sheets) (16 packages of 150 sheets)	
2009828-061	ECG paper, letter, red grid with blank data area, 10 year archival life, z-fold	Case (2400 sheets) (16 packages of 150 sheets)	
22616605	ECG paper, A4, red grid, 8 year archival life, z-fold	Case (1500 sheets) (10 packages of 150 sheets)	
22616611	ECG paper, A4, red grid, 25 year archival life, z-fold		

MAC 2000

Part number	Product description	Case qty	Visual aid
ECG trunk cables			
22341808V	Cable, DCAR, 12-lead for MAC systems, 210 cm	Each	
ECG trunk cables with integrated leadwires			
2029890-001V	Cable trunk MAC500/1200 10-lead	Each	
ECG leadwire sets			
2089797-004	Leadwire set, DCAR universal, mixed length for 12-lead. Set of 10 leads (6 leads 66 cm, 2 leads 90 cm, 2 leads 100 cm long)	Each (set of 10)	
38401816	Leadwire, DCAR, banana-ended, 12-lead Set of 10 leads (4 leads with length of 1295 mm, 6 leads with length of 737 mm)	Each (set of 10)	
Banana connectors for use with universal leadwires			
900178-101	Connector, DCAR, banana for 12-lead	Each (set of 10)	
Universal adapter clips for use with banana-ends			
2056813-010	Adapter clip, DCAR, Sure-Lock clip for 12-lead	Package (10 each)	
2066867-010	Adapter clip, DCAR, wide mouth clip for 12-lead	Package (10 each)	
Mactrode alligator clips for use with universal leadwires			
900178-102	Adapter clip, DCAR, Mactrode alligator clip for 12-lead	Each (set of 10)	
Grabber adapter clips for use with universal leadwires			
900178-103	Connector, DCAR, grabber for 12-lead	Each (set of 10)	
ECG accessories			
401089-003	Accessory, DCAR, leadwire separators	Each (set of 6)	
ECG paper			
2009828-060	ECG paper, letter, red grid, 10 year archival life, z-fold	Case (2400 sheets) (16 packages of 150 sheets)	
2009828-061	ECG paper, letter, red grid with blank data area, 10 year archival life, z-fold	Case (2400 sheets) (16 packages of 150 sheets)	
22616605	ECG paper, A4, red grid, 8 year archival life, z-fold	Case (1500 sheets) (10 packages of 150 sheets)	
22616611	ECG paper, A4, red grid, 25 year archival life, z-fold		

MAC 3500

Part number	Product description	Case qty	Visual aid
ECG trunk cables			
22341808V	Cable, DCAR, 12-lead for MAC systems, 210 cm	Each	
ECG leadwire sets			
2089797-004	Leadwire set, DCAR universal, mixed length for 12-lead. Set of 10 leads (6 leads 66 cm, 2 leads 90 cm, 2 leads 100 cm long)	Each (set of 10)	
38401816	Leadwire, DCAR, banana-ended, 12-lead Set of 10 leads (4 leads with length of 1295 mm, 6 leads with length of 737 mm)	Each (set of 10)	
Banana connectors for use with universal leadwires			
900178-101	Connector, DCAR, banana for 12-lead	Each (set of 10)	
Universal adapter clips for use with banana-ends			
2056813-010	Adapter clip, DCAR, Sure-Lock clip for 12-lead	Package (10 each)	
2066867-010	Adapter clip, DCAR, wide mouth clip for 12-lead	Package (10 each)	
Mactrode alligator clips for use with universal leadwires			
900178-102	Adapter clip, DCAR, Mactrode alligator clip for 12-lead	Each (set of 10)	
Grabber adapter clips for use with universal leadwires			
900178-103	Connector, DCAR, grabber for 12-lead	Each (set of 10)	
ECG accessories			
401089-003	Accessory, DCAR, leadwire separators	Each (set of 6)	
ECG paper			
2009828-021	ECG paper, A4, red grid, 10 year archival life, z-fold	Case (2400 sheets) (8 packages of 300 sheets))	
2009828-020	ECG paper, letter, red grid, 10 year archival life, z-fold		
2009828-034	ECG paper, letter, red grid, 25 year archival life, z-fold		
2009828-035	ECG paper, A4, red grid, 25 year archival life, z-fold		
2009828-060	ECG paper, letter, red grid, 10 year archival life, z-fold	Case (2400 sheets) (16 packages of 150 sheets))	
2009828-061	ECG paper, letter, red grid with blank data area, 10 year archival life, z-fold	Case (2400 sheets) (16 packages of 150 sheets))	
2009828-024	ECG paper, letter, red grid with blank data area, 10 year archival life, z-fold	Case (2400 sheets) (8 packages of 300 sheets))	

MAC 5500 HD

Part number	Product description	Case qty	Visual aid
ECG trunk cables			
2016560-001	Cable, DCAR, CAM-14 cable for use with MAC5500 HD, 1.3 m	Each	
2016560-003	Cable, DCAR, CAM-14 cable for use with MAC5500 HD, 4.6 m	Each	
ECG leadwire sets			
2089797-004	Leadwire set, DCAR universal, mixed length for 12-lead. Set of 10 leads (6 leads 66 cm, 2 leads 90 cm, 2 leads 100 cm long)	Each (set of 10 mixed lengths)	
2089797-003	Leadwire set, DCAR universal, mixed length for 15-lead. Set of 14 leads (10 leads 66 cm, 2 leads 90 cm, 2 leads 100 cm long)	Each (set of 14 mixed lengths)	
38401816	Leadwire, DCAR, banana-ended, 12-lead. Set of 10 leads (4 leads 1295 mm long, 6 leads 737 mm long)	Each (set of 10 mixed lengths)	
2016033-002	Leadwire, DCAR, banana-ended, for 15-lead, mixed lengths. Set of 14 leads (2 leads 660 mm long, 2 leads 914 mm long, 10 leads 1016 mm long)	Each (set of 14 mixed lengths)	
Individual ECG leadwires			
2089799-004	Leadwire, DCAR universal, C3R, 66 cm	Each	
2089799-005	Leadwire, DCAR universal, C4R, 66 cm	Each	
2089799-006	Leadwire, DCAR universal, C7, 66 cm	Each	
Banana connectors for use with universal leadwires			
900178-101	Connector, DCAR, banana for 12-lead	Each (set of 10)	
900178-201	Connector, DCAR, banana for 15-lead	Each (set of 14)	
Universal adapter clips for use with banana-ends			
2056813-010	Adapter clip, DCAR, Sure-Lock clip for 12-lead	Package (10 each)	
2056813-014	Adapter clip, DCAR, Sure-Lock clip for 15-lead	Package (14 each)	
2066867-010	Adapter clip, DCAR, wide mouth clip for 12-lead	Package (10 each)	
2066867-014	Adapter clip, DCAR, wide mouth clip for 15-lead	Package (14 each)	
Mactrode alligator clips for use with banana-ends			
900178-102	Adapter clip, DCAR, Mactrode alligator clip for 12-lead	Each (set of 10)	
900178-202	Adapter clip, DCAR, Mactrode alligator clip for 15-lead	Each (set of 14)	
Grabber adapter clips for use with banana-ends			
900178-103	Connector, DCAR, grabber for 12-lead	Each (set of 10)	
900178-203	Connector, DCAR, grabber for 15-lead	Each (set of 14)	

Part number	Product description	Case qty	Visual aid
ECG accessories			
2001926-001	Accessory, DCAR, leadwire blank plug for CAM-14	Pack (set of 4)	
401089-003	Accessory, DCAR, leadwire separators, set of 6	Each (set of 6)	
ECG paper			
2009828-021	ECG paper, A4, red grid, 10 year archival life, z-fold	Case (2400 sheets) ((8 packages of 300 sheets))	
2009828-020	ECG paper, letter, red grid, 10 year archival life, z-fold	Case (2400 sheets) ((8 packages of 300 sheets))	
2009828-034	ECG paper, letter, red grid, 25 year archival life, z-fold	Case (2400 sheets) ((8 packages of 300 sheets))	
2009828-035	ECG paper, A4, red grid, 25 year archival life, z-fold	Case (2400 sheets) ((8 packages of 300 sheets))	
2009828-024	ECG paper, letter, red grid with blank data area, 10 year archival life, z-fold	Case (2400 sheets) ((8 packages of 300 sheets))	
2009828-061	ECG paper, letter, red grid with blank data area, 10 year archival life, z-fold	Case (2400 sheets) ((16 packages of 150 sheets))	

SEER Light

Part number	Product description	Case qty	Visual aid
ECG leadwire sets			
2008594-003	Leadwire, DCAR Holter, SEER Light, 2CH (5 leadwires), 140 cm	Each	
2008594-004	Leadwire, DCAR Holter, SEER Light, 3CH (7 leadwires), 140 cm	Each	
Accessories			
2008596-001	Carry Case, Belt, Strap, SEER Light	Each	

SEER 1000

Part number	Product description	Case qty
ECG leadwire sets		
2067634-050	Leadwire, DCAR Holter, SEER 1000, 3-lead 3CH, 75 cm	Each
2067634-051	Leadwire, DCAR Holter, SEER 1000, 3-lead 3CH, 105 cm	Each
2067634-013	Leadwire, DCAR Holter, SEER 1000, 5-lead 2CH, 75 cm	Each
2067634-015	Leadwire, DCAR Holter, SEER 1000, 5-lead 2CH, 105 cm	Each
2067634-014	Leadwire, DCAR Holter, SEER 1000, 7-lead 3CH, 75 cm	Each
2067634-016	Leadwire, DCAR Holter, SEER 1000, 7-lead 3CH, 105 cm	Each

KISS system

Part number	Product description	Case qty
ECG leadwires for KISS system		
38401590	KISS leadwire set, 10/set, mixed	1 each (set of 10)
30344365	Leadwire, KISS, 50" (127 cm), C1	Each
30344366	Leadwire, KISS, 50" (127 cm), C2	Each
30344367	Leadwire, KISS, 50" (127 cm), C3	Each
30344368	Leadwire, KISS, 50" (127 cm), C4	Each
30344369	Leadwire, KISS, 50" (127 cm), C5	Each
30344370	Leadwire, KISS, 50" (127 cm), C6	Each
30344371	Leadwire, KISS, 57" (147 cm), F	Each
30344372	Leadwire, KISS, 57" (147 cm), L	Each
30344373	Leadwire, KISS, 57" (147 cm), R	Each
30344374	Leadwire, KISS, 57" (147 cm), N	Each
30344375	Leadwire, KISS, 50" (127 cm), NST	Each
30344376	Leadwire, KISS, 50" (127 cm), NAX	Each
KISS leadwires accessories		
21732801	Electrode, KISS suction, ES600 reusable	Each
93204775	Electrode, KISS filter disc	PAK (100 each)
43252125	Electrode, DCAR reusable suction cup	Each
30344421	Adapter clip for adhesive electrodes	Each
Accessories		
21708318	Cream, electrode, 250 mL	Each
21730702	Spray bottle, electrode 250 mL	Case (10 each)
21730705	Spray, electrode refill 2 L	Each
2034731-002	Electrode gel 250 mL 4 bottles/case	Case (4 each)
2034741-002	Epicont bottle, 250 mL	Case (4 each)



21732801, electrode, KISS, suction



93204775, electrode, KISS, filter disc



43252125, electrode, DCAR reusable suction cup

CM 3000 SM, SM/A, Z, ZSM; CM 3000-12(BT); SEER 12

Part number	Product description	Use with	Compatibility	Case qty
5-lead combined cable & leadwire				
2014606-040	GETEMED AMBUL-EKG PAT-CABLE 1PC DIGI 5P (Patient cable with leadwires, 5 leads, 180 cm)	N/A	CM 3000 SM, SM/A, Z, ZSM	1 Each
2014606-070	GETEMED AMB-EKG PAT-CABLE 1PC 5P 0.8 m (Patient cable with leadwires, 5 leads, 80 cm)	N/A	CM 3000 SM, SM/A, Z, ZSM	1 Each
2014606-185f	GETEMED CABLE 5-LEAD ONE-PART 1.2 m (Patient cable with leadwires, 5 lead, 120 cm)	N/A	CM 3000 SM, SM/A, Z, ZSM	1 Each
5-lead separate cable & leadwire				
2014606-071	GETEMED AMB-EKG-PAT-CABLE SW DIG 5P 0.5 m (Trunk cable adapter for 5 leadwires, 50 cm)	2014606-038 or 2014606-072	CM 3000 SM, SM/A, Z, ZSM	1 Each
2014606-043	GETEMED AMB-EKG-PAT-CABLE TCH DIG 5P (Patient trunk cable adapter for 5 leadwires, 120 cm)	2014606-038 or 2014606-072	CM 3000 SM, SM/A, Z, ZSM	1 Each
2014606-072	GETEMED PRESS-STUD ADP-SET 5PC DIGI 0.3 m (Leadwire set, 5 leads, 30 cm)	2014606-043 or 2014606-071	CM 3000 SM, SM/A, Z, ZSM	1 Each
2014606-038	GETEMED PRESS-STUD ADAPT-SET 5PC DIGI (Leadwire set, 5-leads, 60 cm)	2014606-043 or 2014606-071	CM 3000 SM, SM/A, Z, ZSM	1 Each
7-lead combined cable & leadwire				
2014606-073	GETEMED AMB-EKG PAT-CABLE 1PC 7P 0.8 m (Patient cable with leadwires, 7 leads, 80 cm)	N/A	CM 3000 SM, SM/A, Z, ZSM	1 Each
7-lead separate cable & leadwire				
2014606-044	GETEMED AMB-EKG-PAT-CABLE TCH DIG 7P (Trunk Cable adapter for 7 leadwires, 150 cm)	2014606-075	CM 3000 SM, SM/A, Z, ZSM	1 Each
2014606-075	GETEMED PRESS-STUD ADP-SET 7PC DIGI 0.3 m (Leadwire set, 7 leads, 30 cm)	2014606-044	CM 3000 SM, SM/A, Z, ZSM	1 Each
12-lead for CM 3000-12(BT) and SEER 12				
2014606-635	GETEMED CABLE PTNT 12-LEAD ROUND IEC (Patient cable with leadwires, 10 leads, 130 cm, IEC)	N/A	CM 3000-12(BT), SEER 12	1 Each

Electrodes

Part number	Product description	Environment	Case qty	Visual aid
9623-810P	Silver Mactrode Plus	Resting	10/card, 100/pouch, 1000/CS	
9623-814P	Silver Mactrode Plus	Resting	14/card, 140/pouch, 1400/CS	
402207-810	BabyMAC, infant	Resting	10/card, 50/pouch, 500/CS	
2014778-001	SilverTrace, wet gel, banana connection	Resting	30/pouch, 300/CS	
2014775-001	SilverTrace, foam, rectangle, wet gel	Resting	30/pouch, 300/CS	
2014780-001	SilverTrace, foam, round, wet gel	Stress, holter, telemetry	30/pouch, 300/CS	
900703-805	SilverTrace, foam, oval, wet gel	Stress, holter, telemetry	5/card, 600/CS	
2014783-001	Silvertrace, vinyl/clear tape	Monitoring, holter	30/pouch, 300/CS	

For additional electrodes, please contact your Vyaire representative



GLOBAL HEADQUARTERS

Vyaire Medical, Inc.
26125 N. Riverwoods Blvd.
Mettawa, IL 60045
USA



See product labeling
and dedicated product
brochure



Vyaire Medical Pty Ltd
Level 5, 7 Eden Park Drive
Macquarie Park, NSW, 2113
Australia

vyaire.com

For Europe, Middle East, Africa, Asia, Australia, New Zealand distribution only.

© 2018 Vyaire Medical, Inc. or one of its affiliates. All rights reserved. Vyaire and the Vyaire Medical logo are trademarks or registered trademarks of Vyaire Medical, Inc or one of its affiliates. Medical Devices Class I, Class I Measuring according to Medical Devices Directive 93/42/EEC. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-INTL-1900001